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## Inhalt

<b>Vorwort</b> .....	5
<b>Studien</b> .....	
Übersicht.....	6
Is there evidence that bracing could provide adequate stability for a 47-year-old man with a deficient anterior cruciate ligament to resume downhill skiing? <i>Phys Ther, 2004</i> .....	8
Effect of bracing on the prevention of anterior knee pain— a prospective randomized study, <i>Knee Surg Sports Traumatol Arthrosc, 2004</i> .....	10
Unload it: the key to the treatment of knee osteoarthritis. <i>Knee Surg Sports Traumatol Arthrosc, 2011</i> .....	11
Functional knee brace use effect on peak vertical ground reaction forces during drop jump landing. <i>Knee Surg Sports Traumatol Arthrosc, 2012</i> .....	12
The effect of prophylactic knee bracing on performance: balance, proprioception, coordination, and muscular power. <i>Knee Surg Sports Traumatol Arthrosc, 2011</i> .....	13
<b>Zusammenfassung</b> .....	14
<b>Quellen</b> .....	16

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Investigators / Source / Year of Publication	Title	Design	Purpose	Evidence Informed Practice
<p>Wetherbee E, Pellecchia GL / Phys Ther / 2004</p>	<p>Is there evidence that bracing could provide adequate stability for a 47-year-old man with a deficient anterior cruciate ligament to resume downhill skiing?</p>	<p>Case generated literature review</p>	<p>To search the literature for evidence that bracing improves knee stability. To guide clinical decision making in the choice of a knee brace for a specific patient at hand.</p>	<p>Suggestions made to the patient were to attain an off-the-shelf FKB (functional knee brace) and to return to skiing conservatively (avoiding ice and moguls). 4 selected studies provided answers to the case at hand: FKB wear improves knee stability and decreases the incidence of subsequent knee injury. Additional patient information consisted of: Explanations that bracing does not replace the need to continue an exercise regimen aimed at maximizing muscle performance during functional activities. Caution against developing a false sense of security while wearing the brace. Caution avoiding activities that might put him at risk for further injury. (Details to conjunct studies see page 8)</p>
<p>Van Tiggelen D et al / Knee Surg Sports Traumatol Arthrosc / 2004</p>	<p>Effect of bracing on the prevention of anterior knee pain - prospective randomized study</p>	<p>Prospective Randomized Study</p>	<p>To assess the effectiveness of a dynamic patellofemoral brace in prevention of anterior knee pain syndrome (AKPS).</p>	<p>Use of a dynamic patellofemoral brace is an effective way to prevent the development of anterior knee pain (AKPS) in persons undergoing a strenuous training program: When braced during 6 weeks of basic military training (BMT) 10 of 54 recruits (18,5%) reported anterior knee pain compared to 42 of 113 (37%) in the control group. 1 braced compared to 4 unbraced out of 5 recruits sustained a stress fracture. 1 braced compared to 7 unbraced out of 8 recruits sustained another traumatic lesion. The results of this study support the use of prophylactic patellofemoral bracing, since prevention of overuse injuries is in many cases a primary concern. The mechanism behind this prophylactic effect of bracing on patellofemoral joint remains uncertain: Bracing may alter proprioceptive input and muscular recruitment positively.</p>
<p>Waller C et al / Knee Surg Sports Traumatol Arthrosc / 2011</p>	<p>Unload it: the key to the treatment of knee osteoarthritis</p>	<p>Selective literature review: Research of trends in the treatment of osteoarthritis (OA)</p>	<p>To consolidate and underscore the most salient thinking developed on the subject of aberrant biomechanics as the cause of OA.</p>	<p>Aberrant biomechanics are the primary driver in the progression of osteoarthritis (OA) of the knee. Reestablishing a favorable local mechanical environment may not only delay the requirement for an invasive joint reconstruction procedure but obviate the need entirely. Unloading or offloading the joint will reverse the structural damage: Clinically, patients are likely to experience symptomatic relief even in the short term. Long term, competent hyaline cartilage may re-form if the load distribution is durably maintained. Selected studies suggest weight loss, wedged insoles and knee braces as minimal-invasive modalities to effectively correct biomechanics of joints and target the progression of OA. HTO (high tibial osteotomy) and joint distraction surgery are invasive procedures. Evidence is lacking that HTO is more effective than conservative care. Valgus unloader knee braces demonstrated significantly better pain relief and functional improvement compared to neoprene sleeves or no brace after 6 month of follow up.</p>
<p>Rishiraj N et al / Knee Surg Sports Traumatol Arthrosc / 2012</p>	<p>Functional knee brace use effect on peak vertical ground reaction forces during drop jump landing</p>	<p>Field Trial: Crossover study</p>	<p>To investigate the landing strategies used by non-injured athletes while wearing functional knee braces (FKB, BR condition) during a drop jump task compared with non-injured, nonbraced (NBR condition) subjects. To ascertain whether accommodation to a FKB was possible by noninjured braced subjects.</p>	<p>FKBs offer protection to an ACL-deficient knee. FKBs address concerns about early muscle fatigue, energy expenditure, heart rate and decrease in performance level. The significantly lower group mean peak vertical ground reaction forces (PVGRF) while using a FKB could keep traumatic forces from reaching the ACL until the active neuromuscular restraints are activated to provide protection to knee joint ligaments. Accommodation to FKB is possible after approximately 14.0 hours of brace use.</p>
<p>Baltaci G et al / Knee Surg Sports Traumatol Arthrosc / 2011</p>	<p>The effect of prophylactic knee bracing on performance: balance, proprioception, coordination, and muscular power</p>	<p>Clinical Test: Crossover study</p>	<p>To determine which one of five selected prophylactic knee braces is most effective on functional performance.</p>	<p>Prophylactic knee braces can be used for both healthy subjects and athletes to enhance proprioception, coordination, maximal force and balance. Knee bracing was shown to provide coordination and maximal force through a restriction of range of motion and an enhancement of the proprioception of the knee. Prophylactic knee braces may inhibit specific performance parameters in subjects unaccustomed to wearing braces, which may result in subsequent decrements in athletic performance. It is important to choose the best proper brace for the individual according to his/her fitness level and functional goals. Out of 5 braces 2 produced the best results on performance parameters: The hinged H buttress for support of the knee brace showed the best effect on dynamic balance tests. The Drytex economy hinged knee brace produced the best results in proprioception and maximal force tests.</p>

## Is there evidence that bracing could provide adequate stability for a 47-year-old man with a deficient anterior cruciate ligament to resume downhill skiing? *Phys Ther, 2004*

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### Electromyographic activity in expert downhill skiers using functional knee braces after anterior cruciate ligament injuries.

*Am J Sports Med, 1997*

Purpose: 6 expert downhill skiers who had sustained ACL injuries with different degrees of knee instability were tested to measure possible changes in electromyographic (EMG) activity during downhill skiing without and with a custom-fit brace applied to the injured knee:

- *With the brace, the EMG activity increased in midphase and the peak activity occurred closer to knee flexion in midphase.*
- *With the brace on the injured leg, the uninjured knee was influenced and a decrease in EMG activity was seen in midphase.*
- *Overall, brace-wearers demonstrated beneficial motor control patterns in response to afferent input from knee proprioceptors and from the brace-skin-bone interface to the central nervous system.*
- *With the FKB, the subjects reported feeling more stable and confident*

### Effect of functional bracing on subsequent knee injury in ACL-deficient professional skiers.

*J Knee Surg, 2003*

Purpose: To evaluate the effect of functional bracing on subsequent knee injury in the ACL in professional skiers (skiing 80 days/season). Deficiency was defined by the Lachman or Pivot-Shift Test and a 5-mm KT-1000 manual maximum difference. From the screening of 9410 subjects 180 ACL-deficient skiers were identified.

- *The skiers that were not braced had a higher proportion of subsequent knee injuries than the braced skiers: 2 of 101 braced skiers and 10 of 79 unbraced skiers sustained subsequent knee injuries.*
- *Logistic regression identified bracing status and KT-1000 manual maximum as significant multivariate risk factor for subsequent knee injury.*

### Functional knee braces and dynamic performance: a review.

*Clin J Sport Med, 1997*

Purpose: To examine the effectiveness of functional knee braces (FKBs) used by patients with ACL-injury during dynamic performance tests. A comprehensive literature review included 12 studies comparing braced and unbraced conditions.

- Bracing was found to be advantageous in 3 of 16 maximal effort situations, disadvantageous in 2 of 16 and of no significant effect in 11 of 16.
- In 5 of the 12 studies reviewed, patients with bracing felt more stable, had fewer episodes of “giving way” and believed that wearing a brace enhanced their performance.
- In some studies, initially low performing patients improved their functional performance by up to 50% when braced.
- Caution is advised in the following: Potential false sense of security and subsequent risk of reinjury, trade-off between stabilization and performance; and factors of injury, bracing, surgical and rehabilitative interventions effecting performance with FKBs.

### “Custom-Fit” versus “Off-the-Shelf” ACL functional braces.

*Am J Knee Surg, 2001*

Purpose: 5 patients with unilateral, chronic ACL tears were each tested without a brace and then in each of the three testing conditions with 2 functional knee braces (custom brace, off-the-shelf brace) and taping.

- *Both FKB braces improved anterior knee stability over knee taping.*
- *From isokinetic testing conclusions were made that neither the custom-fit brace nor the off-the-shelf brace compromised the patients’ muscle force and endurance.*
- *After exercise, brace migration was greater for the custom brace.*

## Effect of bracing on the prevention of anterior knee pain— a prospective randomized study *Knee Surg Sports Traumatol Arthrosc, 2004*

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- Numerous reports exist on therapeutic interventions associated with the use of braces as conservative treatment for AKPS. The function of those braces is to improve the patellar tracking and maintain the patellofemoral alignment by stretching thigh lateral structures.
- Braces may affect the knee region not only by mechanical functions but by other mechanisms such as thermal effects, an increased sensory feedback and an altered circulation.
- Preventative patellofemoral bracing may be viewed as a method to help maintain an "ideal" biomechanical environment in order to avoid irritation of the surrounding tissues.
- In addition to a possible mechanical effect, it has been shown that wearing a knee brace has significant consequences on the proprioceptive capacity.
- A significantly smaller IEMG signal of the VMO and VL during braced isokinetic activity compared to a non-braced condition was noted. The combination of the potential anchoring effect of an external device and its anatomical stabilization may reduce the need for muscle activity. This could influence the tracking of the patella and consequently the appearance of AKPS in subjects with poor patellar tracking.
- Bracing may alter the recruitment patterns of surrounding muscles. Isokinetic performances improved in patients with patellofemoral pain when a brace was used.
- It is probable that a combination of the mechanic mechanisms, the proprioceptive properties and the altered muscular recruitment pattern is responsible for the decreased incidence of AKPS when performing a strenuous training program.

## Unload it: the key to the treatment of knee osteoarthritis. *Knee Surg Sports Traumatol Arthrosc, 2011*

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### Gewichtsverlust

- Diet & exercise with individuals with (BMI) of 34 kg/m<sup>2</sup> reported that each 1 kg of weight loss resulted in a 4 kg reduction in knee compressive forces.
- Even a modest 10 lb (4,5 kg) loss of body weight can lower the cumulative compressive forces at the knee by almost 50,000 lb (22,680 kg) for each mile walked.
- Body weight loss of 11 lb (5 kg) lowers the risk of developing OA by over 50%.

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- Wedged insoles in footwear may offer a simple solution to alter mechanical forces on the knee by shifting the tibio-femoral angle.
- Foot orthoses provide mild symptomatic relief with respect to knee pain and stiffness.
- Several RCTs failed to demonstrate durable symptom relief.
- An estimated 5-10 hour daily use of specialized footwear is required.

### Knie Bandagen (knee unloader braces)

- *Valgus unloader braces in various knee OA demonstrated significantly better pain relief and functional improvement compared to neoprene sleeves or no brace after 6 months of follow-up.*

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HTO (high tibial osteotomy); knee joint distraction

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## Functional knee brace use effect on peak vertical ground reaction forces during drop jump landing.

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## The effect of prophylactic knee bracing on performance: balance, proprioception, coordination and muscular power.

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- University of Hartford, Conn. / University of Connecticut (USA)
- Departments of Orthopaedics and Sports Orthopaedic Surgery, Stockholm (Sweden) / School of Human Kinetics and Department of Anatomy and Neurobiology, University of Ottawa (Canada) / Department of Orthopaedic Surgery, Alexandra Hospital, Singapore (Singapore)
- Centre for Physical Medicine and Rehabilitation, Queen Astrid Military Hospital, Brussels / Department of Rehabilitation Sciences and Physiotherapy, Ghent University / Medical Department, Royal Military Academy, Brussels/ Department of Orthopaedics and Rehabilitation, Ghent University ( Belgium)
- St. Vincent ´s Clinic, Darlinghurst (Australia) / Brisbane Orthopedic and Sports Medicine Center, Brisbane (Australia) /Office San Francisco (USA), Department of Orthopaedic Surgery, Harrogate (UK)
- University of British Columbia, Vancouver (Canada) / University of Otago, Dunedin (New Zealand)
- Department of Physiotherapy and Rehabilitation, Hacettepe University, Ankara (Turkey)

## Associations involved

- NATA - National Athletic Trainers' Association (USA)
- APTA - American Physical Therapy Association (USA)
- ESSKA - European Society of Sports Traumatology Knee Surgery and Arthroscopy

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- Physical Therapy- Journal of the American Physical Therapy Association (USA)
- Journal Knee Surgery (USA)
- Clinical Journal of Sports Medicine (USA; Canada; Australia; Asia)
- American Journal of Sports Medicine via Sage Publications, Inc. (USA)
- Journal of Knee Surgery via Slack Publications, Inc. (Europe)
- American Journal of Knee Surgery (USA)
- Knee Surgery Sports Traumatology Arthroscopy -Official Journal of the European Society of Sports Traumatology Knee Surgery and Arthroscopy (ESSKA)

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Subjective Measures	Objective Measures
<ul style="list-style-type: none"> <li>• Pain               <ul style="list-style-type: none"> <li>- yes/no</li> <li>- continuous data</li> <li>- VAS</li> </ul> </li> <li>• Subjective instability (e.g. 'giving way'; yes/no)</li> <li>• Knee functional mobility</li> <li>• Knee AROM (continuous patient reports and behavior)</li> <li>• Patient satisfaction (interval, continuous or dichotomous data)</li> </ul>	<ul style="list-style-type: none"> <li>• Return to sports at the previous performed level (yes/no; time to achieve)</li> <li>• Return to work (yes/no; time to achieve)</li> <li>• Objective instability               <ul style="list-style-type: none"> <li>- Knee: Lachman, Pivot Shift Test, medial/ lateral ligament tests, meniscal grind, etc.</li> </ul> </li> <li>• Swelling (yes/no; grade of edema)</li> <li>• Recurrent injury (yes/no)</li> <li>• Knee mobility (graded manual testing)</li> <li>• Knee ROM (continuous data)</li> <li>• Complications:               <ul style="list-style-type: none"> <li>- sensory deficit</li> <li>- infection</li> <li>- arthrosis /osteoarthritis</li> <li>- allergic reaction</li> <li>- stiffness</li> <li>- muscle atrophy, etc</li> </ul> </li> </ul>

## Quellen

---

### Researchstudien / Reference List:

- a. Wetherbee E, Pellecchia GL. Is there evidence that bracing could provide adequate stability for a 47-year-old man with a deficient anterior cruciate ligament to resume downhill skiing? *Phys Ther.* 2004; 84:274-282.
  - I. Kocher MS et al. Effect of functional bracing on subsequent knee injury in ACL-deficient professional skiers. *J Knee Surg.* 2003 Apr;16(2):87-92.
  - II. Kramer JF et al. Functional knee braces and dynamic performance: a review. *Clin J Sport Med.* 1997;7:32-9.
  - III. Nemeth G et al. Electromyographic activity in expert downhill skiers using functional knee braces after anterior cruciate ligament injuries. *Am J Sports Med.* 1997;25:635-41.
  - IV. Wojtys, EM, Huston LJ. "Custom-Fit" versus "Off-the-Shelf" ACL functional braces. *Am J Knee Surg.* 2001;4:157-62.
- b. Van Tiggelen D et al. Effect of bracing on the prevention of anterior knee pain – a prospective randomized study. *Knee Surgery Sports Traumatol Arthrosc.* 2004;12(5): 434-439.
- c. Waller C et al. Unload it: the key to the treatment of knee osteoarthritis. *Knee Surg Sports Traumatol Arthrosc.* 2011;19:1823-29.
- d. Rishiraj N et al. Functional knee brace use effect on peak vertical ground reaction forces during drop jump landing. *Knee Surg Sports Traumatol Arthrosc.* 2012;20:2405-12.
- e. Baltaci G et al. The effect of prophylactic knee bracing on performance: balance, proprioception, coordination, and muscular power. *Knee Surg Sports Traumatol Arthrosc.* 2011;19:1722-28.

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**BAUERFEIND AG**
Bauerfeind AG