# **EVIDENCE INFORMED PRACTICE**

THE ANKLE JOINT

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SPRUNGGELENK

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# INTRODUCTION: EVIDENCE RELATING TO THE EFFECTIVENESS OF SUPPORTS/ORTHOSES IN PREVENTING INJURIES AND REINJURIES OF THE ANKLE JOINT

## **What are the benefits associated with the use of supports/orthoses?**

Evidence relating to the role of supports/orthoses in preventing injuries and reinjuries of the ankle joint is provided by the specialist literature.

Positive research results and, therefore, evidence relating to the successful use of supports as a preventive measure were identified by using the following search terms: **"ankle AND prevention of injury, reinjury OR prophylaxis AND elastic brace, bracing, functional support, orthoses, tape, taping."**

The matrix on pages 6 and 7 provides an overview of the studies included and a summary of their contents.

The facts and information backed up by scientific research have a significance that should be communicated to specialists and to their clients and patients (evidence based practice/evidence based medicine).

Further details about the studies are provided on the pages following the study matrix.

Notes: The general term "brace" or "bracing" used in the studies applies to the use of both supports (active supports) and orthoses (stabilizing orthoses, functional orthoses and multifunctional orthoses).

The studies included participants with injuries, those with injuries that had healed and healthy individuals.

## **Scientific institutes, sources/specialist literature and outcome measures:**

International research groups at scientific institutes (page 20) set themselves the task of investigating the management of joint injuries in order to draw conclusions about the best possible treatment practice. Medical specialists should be supported through the availability of evidence and guidelines provided by the literature reference list (page 22).

Outcome measures in the form of validated tests, observations and patient reports (page 21) substantiate both the subjective and objective parameters of the findings.

As well as being used for diagnostic purposes, these are also used to evaluate rehabilitative progress and improve communication about functional targets between physicians, therapists, patients and, in many cases, patients' relatives.

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Investigators/Source/ Year of Publication	Title	Design
Seah R, Mani-Babu S/ Br Med Bulletin/2011	<b>Managing ankle sprains in primary care: what is best practice? A systematic review of the last 10 years of evidence</b>	Systematic Review
Verhagen EA et al/ Clin J Sport Med/2000	<b>The Effect of Preventive Measures on the Incidence of Ankle Sprains</b>	Systematic Review
Verhagen EA, Bay K/ Br J Sports Med/2010	<b>Optimising ankle sprain prevention: a critical review and practical appraisal of the literature</b>	Literature Review
Hardy L et al/J Athl Train/2008	<b>Prophylactic Ankle Braces and Star Excursion Balance Measures in Healthy Volunteers</b>	Crossover Study/Laboratory Setting
Olmsted LC et al/ J Athl Train/2004	<b>Prophylactic Ankle Taping and Bracing: A Numbers-Needed-to-Treat and Cost-Benefit Analysis</b>	Analysis/Literature Review
Kerkhoffs GMMJ et al/Cochrane Database of Systematic Review/ 2009	<b>Different functional treatment strategies for acute lateral ankle ligament injuries in adults</b>	Intervention Review
Lardenoye S et al/BMC Musculo- skeletal Disorders/2012	<b>The effect of taping versus semi-rigid bracing on patient outcome and satisfaction in ankle sprains: a prospective, randomized controlled trial</b>	Randomized Controlled Trial
Janssen KW et al/British Journal of Sports Medicine/2014	<b>Bracing superior to neuromuscular training for the prevention of self-reported recurrent ankle sprains: a three-arm randomised controlled trial</b>	Randomized Controlled Trial

## Purpose

## Evidence Informed Practice

Evaluating best practice in managing ankle sprains by comprehensive review of literature (2000–2009).

- For mild to moderate ankle sprains, functional treatment options are statistically better than immobilization.
- Semi-rigid orthoses and pneumatic braces provide beneficial ankle support and may prevent subsequent sprains during high-risk sporting activity.
- Lace-up supports are a more effective functional treatment than elastic bandaging and result in less persistent swelling when compared with semi-rigid ankle braces.

Critical review of the efficacy of preventive measures on the incidence of lateral ankle ligament injuries described in the literature (1980–1998).

- The use of either tape or braces reduces the incidence of ankle sprains.
- The use of tape or braces results in less severe ankle sprains.
- Braces seem to be more effective in preventing ankle sprains than tape.
- It is not clear which athletes benefit more from the use of preventive measures: those with or those without previous ankle sprains.

To establish the effect of preventive measures and assert the optimal prevention strategy for acute lateral ligament injury to the ankle.

- Taping, bracing and neuromuscular training are all effective for the prevention of ankle-sprain recurrences. The combination of either taping or bracing one with training is preferable.
- Preventive effects have been reported in a general athletic population.
- Effect of bracing is amplified when athletes were previously injured.

To determine if prophylactic ankle braces affected multidirectional reach distances during a test of dynamic balance.

- In a comparison between the conditions of with no brace, semi-rigid and lace-up ankle brace, the bracing condition had no effect on any Star Excursion Balance Test measures in healthy volunteers.
- Clinicians can feel confident that the prophylactic use of ankle braces in healthy participants does not disrupt lower extremity dynamic balance during a reaching task.

To examine the effectiveness of ankle taping and bracing in reducing ankle sprains by applying a numbers-needed-to-treat (NNT) analysis to previously published studies.

- Greater benefit is achieved in applying prophylactic ankle taping or bracing to athletes with a history of ankle sprain, compared with those without previous sprains.
- When deciding whether athletes should be taped or braced, the increased cost (taping is 3 times more expensive) and the time-consuming nature of ankle taping must be considered.
- Ankle bracing may be a better way to provide the support necessary to prevent ankle sprains.

To assess different functional treatment strategies for acute lateral ankle ligament ruptures in adults.

- In the control of persistent swelling, a lace-up ankle support achieved significantly better short-term-results than a semi-rigid ankle support, elastic bandage and tape.
- Use of a semi-rigid ankle support resulted in a significantly shorter time to return to work and sport when compared with an elastic bandage.
- Fewer patients reported instability when treated with a semi-rigid support than with an elastic bandage. Tape treatment resulted in significantly more complications when compared with an elastic bandage.

To determine the effect of treatment with tape compared to treatment with a brace on patient outcome and satisfaction in ankle sprains.

- The treatment of acute lateral ankle sprain with a semi-rigid brace leads to fewer complications (allergic contact dermatitis, bullae formation, skin pressure abnormalities).
- Posthoc testing revealed that during the 4-week treatment period patient satisfaction was significantly higher in the group treated with a brace at 3 and 5 weeks. Satisfaction with the tape treatment significantly decreased from week 1 till week 5.
- There is no difference regarding functional outcome and pain.

To evaluate the effectiveness of combined bracing and neuromuscular training, or bracing alone, against the use of neuromuscular training on recurrences of ankle sprain after usual care in athletes following an acute injury to the lateral ankle ligaments.

- Bracing is superior to neuromuscular training for the prevention of self-reported ankle sprain recurrences. Bracing is associated with an added 47% reduction in risk of recurrence versus neuromuscular training.
- Although the current clinical guidelines are vague on the prescription of neuromuscular training and bracing, the study results support the prescription of bracing as single secondary preventive measure.
- Bracing was proven effective when prescribed during sports for 12 months. The prescription period of brace-use in athletes needs to be extended instead of being phased out.

# MANAGING ANKLE SPRAINS IN PRIMARY CARE: WHAT IS BEST PRACTICE? A SYSTEMATIC REVIEW OF THE LAST 10 YEARS OF EVIDENCE BR MED BULLETIN, 2011

The purpose of this 30-page research review article was to supply evidence about the management of ankle sprains in the last decade (2000–2009). 33 articles were selected on the basis of specific inclusion criteria. 5 systematic reviews, 2 RCTs and 1 cross-sectional survey addressed the topic of conservative treatment strategies. The review focused on the following questions in relation to ankle sprains:

- What are the **best decision-making criteria** (clinical decision rules) and test instruments?
- What is the level of evidence relating to **conservative treatment strategies**?
- What is the level of evidence relating to non-conservative treatment strategies?
- Are there indications that ankle sprains can be **fundamentally avoided or prevented** by taking specific preventive measures?

The study also revealed **predisposing factors** for ankle sprains:

- previous history of ankle sprains/recurrent sprains resulting in chronic instability problems
- ligament complex hyperlaxity
- poor sensorimotor control
- biomechanical abnormalities

The incidence of upper and lower ankle sprains is stated to be 1 injury per 10,000 people per day. 25% of all sports injuries are lateral ligament complex injuries of the ankle; an exception to this statistic is the US military, where the incidence rate, at 35%, is extremely high.

The following results indicate preferred strategies:

- The Ottawa Ankle Rules (OAR) are a well-validated decision aid for identifying a suspected fracture of the upper ankle.
- For mild to moderate ankle sprains, **functional treatment options were statistically better than immobilization**. Functional interventions included **elastic bandaging, soft casts, taping and orthoses combined with coordination training**.
- Semi-rigid orthoses and U-Schiennen braces promoted excellent **stability in the ankle and, subsequently, prevention of injuries in high-risk sporting activities**.
- **Lace-up supports/orthoses were a more effective functional treatment option than elastic bandaging**. In the short-term, they prevented swelling more effectively than semi-rigid orthoses, elastic bandaging and tape.
- **There was no evidence that high-top shoes have any protective effect against** reinjury to the ankle.

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# THE EFFECT OF PREVENTIVE MEASURES ON THE INCIDENCE OF ANKLE SPRAINS

## CLIN J SPORT MED, 2000

This review provides a critical insight into selected literature (1980–1998), focusing on the efficacy of preventive measures on the incidence of ankle sprains. From 254 studies, 8 studies whose quality benchmarks matched the inclusion criteria were ultimately selected. The prevention strategies of these studies included the following options:

- injury prevention exercise program
- a low and high cut of shoe
- preventive taping combined with a low or high cut of shoe
- taping versus lace-up braces
- braces (ankle stabilizers)
- semi-rigid braces
- braces combined with disk training
- outer shoe brace

The studies included the following activities, evaluated over a period ranging from one season (minimum) to six years (maximum):

- volleyball
- basketball
- american football
- football
- skydiving

Extracted clinical data about measures that result in a significant reduction in the incidence of ankle ligament injuries:

- The preventive use of either tape or braces **reduces the incidence of ankle sprains**. However, compared to tape, **braces have no limiting effect in various exercise periods**.
- When sprains occur, the use of tape or braces helps to **reduce the severity of the ankle injury**.
- **Braces are more effective at preventing ankle sprains** than tape.
  - In a 6-year retrospective study on 297 college football players, **lace-up ankle braces were twice as effective as taping** at preventing sprains: “2.56 sprains versus 4.91 sprains per 1000 injury exposures”.
- Tape and braces can both be applied independently by the user, but
  - **braces are easier for the user to put on**, and
  - **braces are more cost-effective**.
- The **newness of the footwear** seems to be more important with regard to stability than the height of the cut of the shoe.
- **Proprioceptive training reduces the incidence of ankle sprains in people with previous sprains** to the same degree as in people with no previous history of ankle sprains.

It is not clear from this study which athletes benefit more from the use of these preventive measures: those with or those without a previous history of ankle sprains.

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# OPTIMISING ANKLE SPRAIN PREVENTION: A CRITICAL REVIEW AND PRACTICAL APPRAISAL OF THE LITERATURE BR J SPORT MED, 2010

This article is a literature review that includes RCTs (randomized clinical trials), CTs (clinical trials) and TIs (timed interventions) from 1999 to 2009. **24 studies** met the inclusion criteria, with the prevention strategies including the following options:

- **Preventive taping**
  - combined with either a low or a high cut of shoe
  - fibular repositioning taping
- **Bracing**
  - semi-rigid ankle braces
  - semi-rigid “ankle stabilizers”
  - brace and balance-board training
  - brace versus balance-board training versus strength training
  - laced braces versus tape
- **Preventions-Exercise programs**
  - balance training
  - supervised rehab training including posture training
  - unsupervised home exercise program including neuromuscular training
- **Footwear options**
  - low or high cut of shoe
  - padded or non-padded upper

The studies included the following activities over a period of 1 season (minimum) to 3 seasons (maximum):

- volleyball
- basketball
- american football
- football
- handball
- skydiving

In summary, the researchers reached the following conclusions about optimized strategies for preventing ankle injuries:

- The most effective form of prevention for athletes is a **combination of external prophylactic measures (brace or tape) and neuromuscular training**.
- Preventive measures have a significantly beneficial effect, especially for athletes who have sustained a previous injury and who are in the rehabilitation phase (**reinjury during the sub-acute stage**). In comparison, the effects with regard to new ankle injuries are very low.
  - Researchers talk about an **“unpreventable baseline ankle sprain risk”**: in amateur sport, medical care should **particularly be provided for the secondary prevention of ankle sprains**.
- As a general rule and for practical reasons, preference should be given to the **measures preferred by the athletes themselves**. This produces the best results, because the athlete’s compliance is highest for that particular measure. Nevertheless, the literature should be used in order to inform athletes about the best options.

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- Volleyball
- Basketball
- American Football
- Fußball
- Handball
- Fallschirmspringen

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# PROPHYLACTIC ANKLE BRACES AND STAR EXCURSION BALANCE MEASURES IN HEALTHY VOLUNTEERS

J ATHL TRAIN, 2008

This laboratory study, with a crossover design, focused on two issues and measured the test results with SEBT (Star Excursion Balance Test Performance).

1. Do braces that are worn to prevent an ankle injury influence the user's functional performance capacity?
2. Are there differences between 2 different types of brace – semi-rigid versus lace-up – in terms of dynamic balance?

Data concerning multidirectional reach distances during a dynamic 1-leg test was collected from each of the **36 healthy volunteers (18 male students/18 female students)** with three different brace conditions:

- unbraced stance leg
- stance leg with a semi-rigid brace
  - T2 Active Ankle; Active Ankle Systems, (Jeffersonville, IN)
- stance leg with a lace-up ankle brace
  - Ankle Stabilizer Orthosis; Medical Specialties Inc. (Charlotte, NC)

The volunteers wore their own, usual athletic shoes for all tests. They were permitted to complete a series of practice tests before the series of tests for data collection that comprised nine attempts for each of the three brace conditions (in a randomized order). Reach distance was measured for each of these attempts in eight directions: anterior, anterior-medial, medial, posterior-medial, posterior, posterior-lateral, lateral and anterior-lateral. Each individual completed 27 attempts for the data collection series.

The objective was to move the reach leg as far as possible in a particular direction without losing balance in the process. The stance leg was not permitted to move from its position, the reach leg was not permitted to rest on the floor and the reach leg had to return to the original position. A mistrial was declared if these conditions were not met.

The results gave rise to the following conclusions:

- The wearing of a brace in order to prevent ankle instability or an ankle injury **had no significant effect when functional performance was tested** in an established balance test.
- There were small and insignificant differences between the two types of brace with regard to restriction of movement in a particular direction:
  - The semi-rigid brace restricted frontal ROM to a minimal degree, if at all.
  - The lace-up brace restricted the sagittal ROM to a minimal degree, if at all.
- The preventive use of ankle braces **does not disrupt or impair the individual's dynamic balance or reach.**
- Athletes and coaches can be certain that the preventive use of ankle braces **does not significantly alter the mobility capabilities of the stance leg under dynamic conditions.**

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# PROPHYLACTIC ANKLE TAPING AND BRACING: A NUMBERS-NEEDED-TO-TREAT AND COST-BENEFIT ANALYSIS

## J ATHL TRAIN, 2004

This literature review provides an interesting analysis of the costs and time commitment involved with two preventive measures in relation to ankle sprains. The purpose of the study was to examine the effectiveness of ankle taping and ankle bracing for the prevention of ankle sprains. **Eight studies** were identified; it was ultimately possible to calculate the **NNT (numbers-needed-to-treat) analysis for three studies.**

The three studies all included male test subjects who trained at least six days a week and who could therefore be defined as "athletes":

- Study 1: **College basketball players (2562 player games),**
- Study 2: **Military academy basketball players (n = 1601)** and
- Study 3: **Semi-professional footballers (n = 504).**

The studies compared the following interventions:

- Study 1: Taping versus bracing versus control
- Study 2: Bracing versus control
- Study 3: Sport-Stirrup Orthosis versus control

Each of the three studies distinguished between athletes who had sustained a previous injury (positive history) and those without a previous injury (negative history). The injury rates were documented for both the treatment group and the control group.

Positive History	Treatment Group	Control Group
Study 1	13 sprains/ 794 athlete exposures	24 sprains/ 434 athlete exposures
Study 2	1 sprain/ 87 subjects	6 sprains/ 90 subjects
Study 3	16 sprains/ 127 subjects	32 sprains/ 117 subjects

Negative History	Treatment Group	Control Group
Study 1	4 sprains/ 367 athlete exposures	12 sprains/ 670 athlete exposures
Study 2	10 sprains/ 702 subjects	29 sprains/ 722 subjects
Study 3	32 sprains/ 117 subjects	33 sprains/ 129 subjects

The conclusions of the analysis were as follows:

- For athletes with a history of ankle injury (ankle sprains) in particular, both options, i.e. ankle taping or bracing, are an effective way of preventing reinjury.
- When deciding on which treatment to propose to an athlete (club-level athletes and semi-professional/professional athletes), both the cost and the time involved for a particular treatment are key factors that should be considered:
  - The treatment ratio (bracing versus taping) was 1:78 during a 13-week season. This means that an athlete needed to be re-taped 78 times on average in order to achieve the same sense of security and stability provided by a brace that had to be put on once before participating in sport on each occasion.
  - The costs of taping are 3 times higher than the costs of supplying a brace.

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# DIFFERENT FUNCTIONAL TREATMENT STRATEGIES FOR ACUTE LATERAL ANKLE LIGAMENT INJURIES IN ADULTS

## COCHRANE DATABASE OF SYSTEMATIC REVIEWS, 2009

This 47-page review is based on the principle of early mobilization as a preferable treatment strategy for acute lateral ankle ligament ruptures. It assesses different functional treatment strategies, and analyzes and compares the following options:

- elastic bandage versus Tape
- elastic bandage versus Semi-rigid ankle support
- elastic bandage versus Lace-up ankle support
- tape versus Semi-rigid ankle support
- tape versus Lace-up ankle support
- semi-rigid ankle support versus Lace-up ankle support

A search was conducted to find studies with the following variable definitions of interventions:

- **“Elastic bandage/stocking”**  
including all types that provide support using an elastic sock-like material (i.e. MalleoTrain)
- **“Tape”**  
including all types that provide support using adhesive and elastic athletic tape
- **“Lace-up ankle support”**  
including all types that provide support using a soft canvas-like or nylon material (i.e. Push-brace)
- **“Semi-rigid ankle support”**  
including all types that provide support through a firm thermo-plastic material comprising a stirrup or posterior rigid support (i.e. Aircast Sport-Stirrup)

By applying specific selection criteria, nine randomized clinical trials (RCTs) with a total of 892 participants were ultimately included.

Significant results can be summarized into the following consensus:

- **Lace-up ankle supports control swelling/edema better** than semi-rigid ankle supports, elastic bandages or tape.
- In terms of short-term results, **semi-rigid supports** were shown to be **significantly better** than elastic bandages in the following categories:
  - faster **return to work**
  - faster **return to sports activities**
  - **less feeling of instability and of shorter duration**
- In the “few complications” category, elastic bandages come out on top, but the recovery process can be achieved more quickly with a semi-rigid or lace-up support. Treatments involving taping resulted in significantly more complications (skin irritation etc.)
- In this review, it was not possible to reach a conclusion about which treatment is the most effective, both clinically and in terms of costs.

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# THE EFFECT OF TAPING VERSUS SEMI-RIGID BRACING ON PATIENT OUTCOME AND SATISFACTION IN ANKLE SPRAINS: A PROSPECTIVE, RANDOMIZED CONTROLLED TRIAL

## BMC MUSCULOSKELETAL DISORDERS, 2012

The researchers investigated the treatment of **acute/sub-acute Grade II and III ankle ligament ruptures caused by inversion trauma**. The clinical trial was conducted in a trauma outpatient clinic in the Netherlands. Patients who presented in the outpatient clinic for secondary treatment following an acute lateral ankle injury (predominantly the ATFL) 5-7 days after the trauma were divided into two groups:

- Patients in one group were taped with Coumans-bandage (3-layer rigid tape bandage that had to be reapplied at least once after two weeks or as required);
- The other group was treated with a semi-rigid brace (AirLoc by Bauerfeind, Zeulenroda, Germany).

Patients in both groups were given instruction and supervision in sensorimotor, range of motion and strength training exercises. All patients were given a home exercise program to be performed daily for four weeks.

Between the original injury and the start of the study, the patients were treated in accordance with the standard RICE protocol (rest, ice, compression and elevation), together with a compressive bandage. Pain medication and the use of crutches were not a standard treatment; instead these were optional extra measures.

Double blinding was applied: Treating physicians and patients were allocated to groups in a randomized process and tests (outcome measures) devised without knowledge of the treatment method. The treatment for both groups lasted for four weeks and follow-up tests were performed in weeks 2, 4, 8 and 12 after the start of the study. Patients were excluded from the study if they had a Grade I inversion trauma, a fracture or a previous ankle injury and if their age was under 16 or over 55 years.

Primary evaluations (primary outcome measures) included patient satisfaction (measured by verbal rating scale), skin or pressure complications (such as allergic contact dermatitis, bullae formation, skin pressure abnormalities) and pain (5-point pain scale).

Secondary evaluations (secondary outcome measures) included ankle function in the form of the Karlsson scoring scale and range of movement measurements performed using a goniometer.

**The RCT ultimately included 100 patients with Grade II and III lateral ankle inversion trauma injuries that were treated either with tape or an ankle brace.** The trial began 5-7 days after the initial acute injury and lasted for four weeks with two follow-up tests in weeks 8 and 12. 38% of the patients had sustained these injuries during sport or sport-related activities.

- **Patient satisfaction with the treatment method was significantly higher in the group with the semi-rigid brace** (especially in week 4), whereas patient satisfaction in the tape treatment group decreased steadily from week 1 to week 5.
- **59.1% of the patients in the group treated with tape experienced complications in the form of skin allergies, pressure problems, hygiene issues and instability problems.**
- Patients in both groups received another type of functional (as opposed to immobilizing, passive) treatment method. **Treatment with a functional focus can be rated as positive in the sense that no differences were identified between the two groups regarding the function regained.**
- **Both functional treatment methods** enabled an identical time of the return to sport or to normal work and both the passive and active ROM were similarly good in both groups.
- From an acute medicine perspective, the use of a brace for a duration of four weeks is more expensive than treatment with tape. On average, the total cost for treatment with a brace was \$204 compared to \$167 for taping.

However, the study reached the conclusion that, **although a higher level of comfort can lead to higher treatment costs, there are also fewer possible complications in comparison to taping.** It is further commented in the study that the short duration of the trial does not allow for the equalization of costs in the cost analysis.

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# BRACING SUPERIOR TO NEUROMUSCULAR TRAINING FOR THE PREVENTION OF SELF-REPORTED RECURRENT ANKLE SPRAINS: A THREE-ARM RANDOMISED CONTROLLED TRIAL

## BRITISH JOURNAL OF SPORTS MEDICINE, 2014

This article discusses a randomized controlled trial in which **384 athletes aged between 18 and 70 years, who had sustained an ankle injury**, were divided into three groups:

- NTG: The training group received instruction in a neuromuscular training program (balance board, exercise sheets and an instructional DVD). Participants were expected to complete 3 training sessions each week for eight weeks, with each session lasting a maximum of 30 minutes. The exercises were made slowly and progressively more difficult with an increasing level of loading.
- BG: The brace group was provided with a semi-rigid brace (Aircast A60, Ankle Support, DJO, Europe) and instructed to wear it during all sports activities for twelve months.
- CG: The combi group received both a semi-rigid brace to be worn during all sports activities for eight weeks and a neuromuscular training program.

The participants included in the study had sustained a lateral ankle injury not more than two months prior to the start of the study, were aged between 18 and 70 years and actively participated in sport for at least 1 hour per week. These individuals responded to advertisements and Internet announcements and were assessed for their suitability for inclusion in the study by means of an interview (oral assessment). Web-based questionnaires were completed once a month for twelve months. Recurrent ankle injuries were documented in an injury registry.

The researchers reached the following conclusions:

- **Bracing is superior to neuromuscular training for the prevention of self-reported ankle sprain recurrences.**
- **Bracing** and, ranked in second place, the combination of neuromuscular training and the wearing of a brace **are a superior rehabilitative measure for reducing the incidence** but not the severity of an ankle sprain recurrence.
  - During the 12-month follow-up, 20% of the study participants (n = 69) reported an ankle sprain recurrence: 29 participants in the NTG, 17 in the BG and 23 in the CG.
  - Per 1000 hours of sport, the incidence of ankle sprain recurrence was 2.51 (NTG), versus 1.78 (CG), versus 1.34 (BG).
- **Semi-rigid braces are associated with a 47% reduction in risk.**
- The study results **support prescription of a brace as the single most effective measure for preventing the recurrence of ankle sprains.**

Instead of gradually reducing the duration for which a brace is worn and ultimately phasing out its use altogether, it **is advocated that its use during sports activities should be extended (up to 12 months).**

Additional background information and critical remarks about the study design: The term "athlete" is used to describe test subjects who participate in sport for at least one hour per week. The discussion addresses the possible limitations posed by concomitant treatments ("intervention spill over"), the problem of self-reported information (web-based questionnaires) and variable compliance with the program. Nevertheless, **bracing is praised in this study as the best-accepted intervention.**

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# SUMMARY

Evidence of the role of joint supports such as support bandages, splints and similar aids following injury to the extremities is supplied by contemporary science and research. There has been extensive research into the functional support provided by the use of supports/orthoses in relation to the ankle in particular. Manufacturers of products designed to support the joints are committed to optimizing these support options for joints and injured structures and, therefore, to increasing the potential of functional rehabilitation.

Pain relief (gate control), the intrinsic drainage of fluids associated with inflammation (edema control) and passive and active joint stability are the key principles and objectives of therapeutic interventions involving the use of supports and orthoses.

The identified studies offer an insight into research into the use of supports/orthoses in the literature review, as well as into the latest scientific findings. Conclusions about the beneficial effect on the ankle of supports/orthoses, both in terms of functional stabilization after injuries and reinjuries, and the prevention of such injuries, have been presented and confirmed.

The studies predominantly emphasize that supports/orthoses make a substantial contribution to the patient's/client's subjective sense of wellbeing. They also support them in their everyday activities and in their efforts to re-achieve their functional targets and sporting performance.

- **Fear of reinjury**, which is generally associated with a protective posture, can be reduced through use of a functional support/orthosis. The wearer feels "safer" during movement.
- Functional supports/orthoses allow activities to be performed within a range of motion (ROM) free from pain and they also support **economization of movement processes**.
- Whether in everyday activities, work, recreational sport or high-performance sport, the prevention of injuries is key to achieving **long-term, high-quality performance**. Performance can be improved through continued practice and uninterrupted training.
- Functional supports/orthoses can reduce the risks of a reinjury and **ensure passive and active stability**, as has been proven by both subjective reports and objective functional tests.
- The **pain management** afforded by wearing supports/orthoses plays a significant role in rehabilitation and **boosts the wearer's mobility. The metabolism of the injured structure is therefore actively and passively stimulated**.
- **Supports and orthoses:**
  - **support** fundamental performance of the structures to be rehabilitated, ensuring this is as physiological and as economical as possible
  - **relieve** the injured structure (**principle of unloading/offloading**)
  - simultaneously relieve the healthy distal and proximal joints and chains of muscles
  - protect bilateral physical structures against degenerative forces that occur due to incorrect posture or misdirected motor unit recruitment
  - free the kinematics as a whole

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# SOURCES

## Scientific institutes responsible for producing the studies included

- Imperial College Healthcare NHS Trust, London (UK)
- Department of Public and Occupational Health, EMGO-Institute, VU University Medical Center, Amsterdam (The Netherlands)
- Western Michigan University, Indiana State University (USA); Calgary, Alberta (Canada)
- The Pennsylvania State University – Department of Kinesiology, PA (USA)
- Department of Orthopedic Surgery, Academic Medical Center, Amsterdam, The Netherlands/Department of Public Health and Primary Care, Leiden University Medical Center, Leiden (The Netherlands)
- Maastricht University Medical Center – Division of Trauma Surgery/Department of Epidemiology (The Netherlands)

## Sources/specialist literature of the studies included

- British Medical Bulletin (UK)
- Clinical Journal of Sport Medicine (Canada)
- British Journal of Sports Medicine (UK)
- Journal of Athletic Training (USA)
- Cochrane Database of Systematic Reviews (International Network)
- BioMed Central Musculoskeletal Disorders (Open Access Online Journal)

## Outcome measures in selected studies

Subjective Measures	Objective Measures
<ul style="list-style-type: none"> <li>• Pain               <ul style="list-style-type: none"> <li>– yes/no</li> <li>– continuous data</li> <li>– VAS</li> </ul> </li> <li>• Subjective instability (e.g. 'giving way'; yes/no)</li> <li>• (Longer) use of assistive devices</li> <li>• Ankle functional mobility (abilities of performance, physical activity participation)</li> <li>• Ankle AROM (continuous patient reports and behavior)</li> <li>• Patient satisfaction (interval, continuous or dichotomous data)</li> <li>• Patient reported hygiene</li> <li>• Cessation of treatment/path</li> <li>• Compliance with treatment/path</li> </ul>	<ul style="list-style-type: none"> <li>• Return to sports at the previous performed level (yes/no; time to achieve)</li> <li>• Return to work (yes/no; time to achieve)</li> <li>• Objective instability Ottawa Ankle Rules; anterior drawer test, talar tilt etc.</li> <li>• (Longer) need for assistive devices</li> <li>• Swelling (yes/no; grade of edema)</li> <li>• Recurrent injury (yes/no)</li> <li>• Ankle mobility (graded manual testing)</li> <li>• Ankle ROM (continuous data)</li> <li>• Karlsson Scoring Scale pain, swelling, instability, stiffness, stair climbing, running, work activities, support</li> <li>• Complications:               <ul style="list-style-type: none"> <li>– (chronic) edema</li> <li>– sensory deficit</li> <li>– arthrosis/osteoarthritis</li> <li>– stiffness</li> <li>– motor deficit/muscle atrophy</li> <li>– skin irritations: allergic contact</li> <li>– dermatitis, bullae, skin pressure</li> <li>– abnormalities etc.</li> </ul> </li> </ul>

## Research Studies

- Shea R, Mani-Babu S., Managing ankle sprains in primary care: what is best practice? A systematic review of the last 10 years of evidence. Br Med Bull. 2011;97:105-35.
- EA Verhagen et al. I. The effect of preventive measures on the incidence of ankle sprains. Clin J Sport Med. 2000;10(4):291-6. II. Optimizing ankle sprain prevention: a critical review and practical appraisal of the literature. Br J Sport Med. 2010;44:1082-8.
- Hardy L et al. Prophylactic Ankle Braces and Star Excursion Balance Measures in Healthy Volunteers. J Athl Train. 2008;43(4):347-351.
- Olmsted LC et al. Prophylactic Ankle Taping and Bracing: A Numbers-Needed-to-Treat and Cost-Benefit Analysis. J Athl Train. 2004;39(1):95-100.
- Kerkhoffs GMMJ et al. Different functional treatment strategies for acute lateral ankle ligament injuries in adults. Cochrane Database of Syst Rev. 2009. Art. No.: CD002938. DOI: 10.1002/14651858.D002938.
- Lardenoye et al. The effect of taping versus semi-rigid bracing on patient outcome and satisfaction in ankle sprains: a prospective, randomized controlled trial. BMC Musculoskeletal Disorders. 2012;13:81.
- Janssen KW et al. Bracing superior to neuromuscular training for the prevention of self-reported recurrent ankle sprains: a three-arm randomised controlled trial. Br J Sports Med. Epub ahead of print: [2014 Jan 7] doi:10.1136/bjsports-2013.092947.

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Bauerfeind USA, Inc.