

PRODUCT CATALOG

MEDICAL COMPRESSION STOCKINGS

2025 / 2026

BAUERFEIND COMPRESSION EXPERTISE: FOUR DIMENSIONS OF EXCELLENCE

Bauerfeind's comprehensive compression expertise ensures that patients across the globe benefit from improved health and mobility – staying true to the motto:

**By people for people.
For a better quality of life.**



POWER OF COMPRESSION

Bauerfeind uses its in-depth medical expertise to develop effective compression solutions for treatment.

ENGINEERED COMPRESSION

Through high-precision engineering, Bauerfeind creates innovative compression products built to the most exacting technical standards.

INNOVATION FOR COMPRESSION


As an innovation leader, Bauerfeind establishes new standards in compression technology, delivering advanced solutions for patients.

BAUERFEIND INNOVATION CENTER

Rigorous testing and quality controls in the Bauerfeind Innovation Center guarantee the reliability and effectiveness of each product generation.

SYSTEMATIC APPROACH TO COMPRESSION FOR THE RIGHT CHOICE

The Bauerfeind systematic approach will help you find the right VenoTrain compression stocking for your customers. The right choice depends on the severity of the condition.



¹ Lurie F. et al. The 2020 update of the CEAP classification system and reporting standards. J Vasc Surg Venous Lymphat Disord, 2020;8(3):342-52.

² Rabe, E. et al.: Bonner Venenstudie der DGP, Phlebologie 1 / 2003

BAUERFEIND.

BY PEOPLE FOR PEOPLE.

FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

We have been developing compression products for phlebological and lymphatic conditions for decades, approaching our production and service with untiring dedication and our development with ingenuity.

EFFECTIVE VARIETY

Our medical **VenoTrain compression stockings** cover all uses/indications in modern compression therapy, offering medically proven effectiveness and unique material properties. From lifestyle and sports products to the treatment of lymphedema: our indication-specific solutions are tailored to individual needs and include a wide variety of circular-knit and flat-knit designs that are available as standard or custom-made solutions.

UNCOMPROMISING QUALITY

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements when it comes to quality, effectiveness and wearing comfort. This is confirmed by inspection procedures, for which our products are permitted to bear the seals and marks after successful testing.

All medical VenoTrain compression stockings are CE-certified, thus complying with current applicable statutory regulations. They also bear the RAL quality label as well as the OEKO-TEX® STANDARD 100 quality seal.

From the thread to the completed product, our medical compression stockings are developed and manufactured in Germany.





VenoTrain® COMPRESSION STOCKINGS FOR SYSTEMATIC CARE

A clear system for optimized patient care

Bauerfeind's well-structured VenoTrain portfolio supports medical retailers in providing optimum patient care. The clear segmentation based on indications and individual requirements ensures simple, fast and accurate product selection, offering numerous advantages for medical retailers:

+ Efficient selection

The straightforward system ensures quick selection of the right compression stocking to suit individual needs.

+ Strong selling points

Each product category highlights specific advantages, helping the medical retailer to effectively communicate its benefits.

+ Tailor-made solutions

The VenoTrain range offers a variety of options for all target groups – from standard products to custom-made solutions.

SIMPLE. CLEAR. CONSISTENT.

The clear structure ensures efficient product selection for optimized care.

THE FIVE PILLARS OF THE VenoTrain® PORTFOLIO

VenoTrain Classic

Comprehensive treatment solutions for patients with chronic venous disorders, ranging from mild to severe symptoms (CEAP stages C0 to C5).

VenoTrain Selection

An exclusive collection of knee-high stockings that combine outstanding materials and style to meet individual needs.

VenoTrain Special

Highly specialized solutions for patients with advanced venous disease and special needs.

VenoTrain Basic

Budget-friendly basic product for chronic venous disease.

VenoTrain Active

Support stockings designed for work, leisure and exercise – perfect for outdoor activities and to combat tired legs.

VenoTrain portfolio



CONTENTS

VenoTrain®
**YOUR LEGS.
YOUR FREEDOM.**

Compression as
active and multi-faceted
as life.



**The right treatment
for every indication.**

CLASSIC

- 06 Product Guide
- 08 **VenoTrain® discrétion**
- 10 Fashion Colors 2025/2026
- 12 **VenoTrain® micro**
- 14 Batik and Design Edition
- 16 **VenoTrain® soft**
- 18 **VenoTrain® soft S**
- 20 Compression therapy in
orthopedics
- 22 **VenoTrain® impuls**

SELECTION

- 24 Product Guide
- 26 **VenoTrain® business**
- 28 **VenoTrain® hero**
- 30 **VenoTrain® cocoon**





SPECIAL

- 32 Product Guide
- 34 **VenoTrain® angioflow**
- 38 **VenoTrain® ulcertec**



BASIC

- 40 Product Guide
- 41 **VenoTrain® pure**



ACTIVE

- 42 Product Guide
- 44 **VenoTrain® act**
Cotton
Sheer Elegance
Microfiber
- 46 **Compression Socks**
Performance
Training

ADDITIONAL PRODUCTS

- 48 **Skin care and other additional products**
- 50 **VenoTrain® glider**
- 51 **VenoTrain® glider plus**

BAUERFEIND-SERVICE

- 52 Measurement Technology
- 53 Academy
- 54 Customer Service
- 55 Online shop

PRODUCT SELECTION AND ORDER

- 56 Edging and top bands
- 58 Circular knit product range
- 60 Additional services | Alterations | Repairs
- 62 How to take correct measurements
- 62 Additional Bauerfeind product categories


CLASSIC

VenoTrain[®] CLASSIC PORTFOLIO



Product guide

The VenoTrain Classic portfolio provides **indication-based solutions** for the treatment of patients with chronic venous disorders. This portfolio covers a wide spectrum: it is aimed at patients with **mild to severe symptoms**, and it includes highly elastic as well as strong compression knits, offering stockings with low or high working pressures. It therefore guarantees **ideal support for CEAP stages C0 to C5** of venous disorders.












Choosing the right Classic product

1. indication and diagnosis provided by a physician
(product specified on prescription)
2. Version and compression class
(product specified on prescription)
3. Patient condition
4. Condition of the connective tissue

Color overview

VenoTrain Classic portfolio

SIMPLE. CLEAR. CONSISTENT.

STANDARD COLORS		micro	soft	soft S	impuls
Caramel		●	●	●	●
Navy		●	●	●	●
Anthracite		●	●	●	●
Black		●	●	●	●
ADDITIONAL COLORS					
Cream		●	●	●	
Amaretto		●	●	●	
Espresso		●	●	●	
FASHION COLORS					
French Lavender		●	●	●	●
Irish Green		●	●	●	●

SIZING SYSTEM

PERFECT-FIT*

VenoTrain micro
VenoTrain soft / soft S

DID YOU
KNOW?

This icon shows which products are available with the Perfect-Fit sizing system.



		XS for VenoTrain soft S		S		M		L		XL	
		CIRCUMFERENCE		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
		cT		up to 85		up to 90		up to 100		up to 120	
		cH		up to 95		up to 105		up to 120		up to 140	
Length cG	short	Silicone top band, wide**		51–57	59–66	56–63	62–69	60–67	66–74	63–71	71–80
	long	Silicone top band, regular	44–50	45–51	52–59	49–56	55–62	53–60	58–66	55–63	62–71
	Type	Top band, AT, silicone top band	44–50	45–57	52–66	49–63	55–69	53–67	58–74	55–71	62–80
		cF	39–50	40–52	43–55	43–55	46–60	45–59	50–64	49–63	52–68
		cE	32–38	33–40	36–44	35–43	39–49	38–48	42–53	42–53	47–58
		cD	28–31	29–35	32–38	32–38	35–42	34–42	37–45	37–45	41–49
		cC	29–33	31–37	34–40	34–40	37–44	37–45	39–47	40–48	43–51
		cB1	23–26	25–29	27–31	27–32	29–34	29–35	31–37	32–38	34–40
		cB	18–20	20–22		22–25		25–28		28–31	
		cY	27–31	29–33		31–35		34–40		37–44	
		cA	17–22	19–24		22–27		25–31		28–35	

* The VenoTrain impuls is measured according to the highly differentiated V6+ sizing system. The VenoTrain discrétion also has a separate sizing system (page 8).

** Wide silicone top band with silicone dots available for: VenoTrain soft, VenoTrain soft S

Sizing system information given in cm

For shoe sizes smaller than 36 or larger than 46 as well as in cases of particular anatomical features of the foot, the use of custom-made products is recommended.

V6+
SIZING
SYSTEM ON
PAGE 22



STANDARD SIZE CHART: VenoTrain discrétion

LEG CIRCUMFERENCE

Size	1	2	3	4
cG	42–56	49–63	56–70	63–77
cC	28–36	32–40	35–43	38–46
cB	18–21	21–24	24–27	27–31

LENGTH

Size	short	long
AG AT	< 70	> 70
AD	< 38	> 38

Sizing system information given in cm

VenoTrain® discrétion

THE TRANSPARENT STOCKING

The exceptionally fine material of these compression stockings gives them a stylish transparency and a smooth shine.

Available in various designs and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure and everyday life.

COLORS

CLAIR



BEIGE



ESPRESSO



BLACK



JAVA



FEATURES

- highly elastic, making the stockings easy to put on
- shapes the legs without constricting them
- breathable fine knitting for a long-lasting wearing comfort

DETAILS

SIZING SYSTEM VenoTrain discrétion (I–IV)

COMPRESSION CLASSES

- Class 1 (10–15 mmHg)
- Class 2 (15–20 mmHg)

TOE closed, open

MATERIAL CONTENT

- Class 1: 75 % polyamide, 25 % elastane
- Class 2: 70 % polyamide, 30 % elastane

CARE INSTRUCTIONS



● STANDARD

	AD Knee stockings	AG Thigh-high stockings with silicone dot pattern band	AG Thigh-high stockings with silicone dotted band	AT Pantyhose
EDGING				
SILICONE TOP BAND				
Clair	●*	●*	●*	
Beige	●	●	●	●
Espresso	●	●	●	●
Black	●	●	●	●
Java	●	●	●	●

*only available in classe 2 (15–20 mmHg)
available in selected countries

experts.bauerfeind.com/venotrain-discretion

OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
598-0153 HOHENSTEIN HTTI
Tested for harmful substances.
www.oeko-tex.com/standard100



C0

C1

C2

C3

C4

C5

C6

VenoTrain®

FASHION COLORS

Brighten up your life.

FRENCH LAVENDER

French Lavender epitomizes strength and sophistication. This expressive color combines elegance with a touch of extravagance – perfect for those who want to embrace their individuality and passion for life through powerful colors.



IRISH GREEN

Irish Green makes a bold statement full of energy and vitality. This fresh, dynamic color highlights the wearer's uniqueness and adds an unmistakable, positive charm to any outfit – a true reflection of individuality and zest for life.



Fashion Colors turn your legs into an eye-catcher and help you pay attention to your health whilst making a fashion statement.

Fashion Colors
available for
VenoTrain micro
VenoTrain soft
VenoTrain soft S
VenoTrain impuls



Visit our Bauerfeind online shop
for more informations:
[experts.bauerfeind.com/
b2b/fashioncolors](https://experts.bauerfeind.com/b2b/fashioncolors)





VenoTrain® micro

THE FASHIONABLE STOCKING

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of the most popular compression stockings.

The rich variety of designs and colors allows you to accentuate your customers' unique style and create a fashion statement.

FEATURES

- Breathable, skin-friendly and easy-care thanks to over 50% microfiber
- Figure-shaping panty area (pantyhose) and a comfort sole for a defined silhouette and high wearing comfort
- High-quality semi-transparent design for an elegant appearance
- Comfort for pregnant women: also available as pantyhose with a special panty area that expands as the mother's bump grows (ATU)

DETAILS



SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 65 % polyamide, 35 % elastane
- Ccl 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS



STANDARD COLORS



ADDITIONAL COLORS



FASHION COLORS



• STANDARD SIZES ▲ MADE-TO-MEASURE SOLUTIONS

	AD	AG	AT	ATU
TOP BAND EDGING				
Top band / Comfort band	● ▲		● ▲	● ▲
Lace sensitive		● ▲		
Silicone micro dots		● ▲		
Silicone dots	● ▲	● ▲		
Silicone dot pattern		● ▲		

VenoTrain® micro

BATIK AND DESIGN EDITION

The Design Edition and the custom-colored Batik variants of the VenoTrain micro will allow your customers to take their style to a new level.

Whether discreetly patterned with a subtle Tango or Jive pattern*, or wildly textured in a Batik design, combined with navy or black, your customers’ legs will undoubtedly stand out as striking eye-catchers.



COLORS

STANDARD COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CARAMEL				
ANTHRACITE				
NAVY				
BLACK				

ADDITIONAL COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CREAM				
AMARETTO				
ESPRESSO				




FASHION COLORS

	TANGO*	JIVE*	BATIK BLACK	BATIK BLACK
FRENCH LAVENDER				
IRISH GREEN				

FRENCH LAVENDER TANGO

DESIGN EDITION

● STANDARD SIZES ▲ MADE-TO-MEASURE SOLUTIONS

	 AD	 AG	 AT
TOP BAND EDGING			
Top band / Comfort band	● ▲		● ▲
Lace sensitive		● ▲	
Silicone micro dots		● ▲	
Silicone dots**	● ▲	● ▲	
Silicone dot pattern		● ▲	

* The Design Edition is available only with closed toes

** not available as Design Edition



VenoTrain® soft ARM STOCKING

Thanks to its comfortable inner elbow zone, the **VenoTrain® soft arm** ensures greater freedom of movement, wearing comfort, and a reduced feeling of pressure.

Length /G

short 39–44

long 44–50

Type CG

cG →

cE →

cC →

CIRCUM- FERENCE

	S min. – max.	M min. – max.	L min. – max.	XL min. – max.
cG	24–30	27–33	29–36	32–40
cE	22–26	24–29	27–32	30–35
cC	14–15	15–16,5	16,5–18	18–20

STANDARD COLORS

BEIGE

CARAMEL

BLACK

DETAILS

COMPRESSION CLASSES

- Ccl 1 (15–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOP BAND EDGING

- Top band
- Silicone dots

More information about the VenoTrain soft arm can be found in the Lymphology Product



VenoTrain® soft

THE ALL-ROUNDER

VenoTrain soft is always the perfect choice: this versatile stocking is designed to withstand even the toughest conditions while remaining smooth and effective.

The all-rounder is available in a variety of designs, offers proven medical effectiveness and provides a soothing massage effect.

FEATURES

- Unisex design
- Gentle on the skin thanks to comfortable materials
- Flexible options – e.g. a single standard stocking in the color caramel
- Comfortable for pregnant women thanks to a special panty area that expands as the mother's bump grows (ATU)

• STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

STANDARD COLORS



ADDITIONAL COLORS



FASHION COLORS



DETAILS

SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open, no foot (made-to-measure), stump (made-to-measure)

MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS



TOP BAND EDGING	AD	AG	AGHB	AT	ATU	ATH	ATX	ATE
Top band / Comfort band	● ▲			● ▲	▲	▲	▲	▲
Lace sensitive		● ▲						
Silicone micro dots		● ▲						
Silicone dots	● ▲	● ▲						
Silicone dot pattern		● ▲						
Garter belt			● ▲					

experts.bauerfeind.com/venotrain-soft



C0

C1

C2

C3

C4

C5

C6



VenoTrain® soft S

THE MULTI-TALENTED STOCKING

Thanks to its spacious Easy heel and reduced knitted fabric in the instep area, our VenoTrain soft S offers particular foot comfort. The anatomically shaped 90° heel and the robust knitted fabric facilitate donning and doffing the product.

This makes the VenoTrain soft S ideal for orthopedic applications, including post-traumatic or post-operative edema management* as well as for people who spend long hours sitting or standing.

STANDARD COLORS



ADDITIONAL COLORS



FASHION COLORS



FEATURES



- Easy heel for easy donning
- High level of wearing comfort thanks to crease-free instep area
- Vitalizing effect thanks to breathable materials

● STANDARD SIZES ▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG
Top band	● ▲	
Lace sensitive		● ▲
Silicone micro dots		● ▲
Silicone dots	● ▲	● ▲
Silicone dot pattern		● ▲

DETAILS



SIZING SYSTEM PERFECT-FIT (XS–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS



* For a complete list of all indications, please see the VenoTrain instructions for use.

experts.bauerfeind.com/venotrain-soft-s



C0

C1

C2

C3

C4

C5

C6

BAUERFEIND – A CONVINCING TEAM

COMPRESSION PRODUCTS FOR ORTHO- PEDIC CARE WITH VenoTrain® soft S

HARNESSING
SYNERGIES



POST-OPERATIVE AND POST-TRAUMATIC KNEE TREATMENT

For targeted treatment following
orthopedic procedures or knee
injuries

RECOM-
MENDED BY
EXPERTS



“Following cruciate ligament surgery, effective stabilization and edema reduction in particular will ensure rapid healing. Our preferred products are orthoses that remain securely in position and medical compression stockings where donning of the product is easy.”

Prof. Dr. med. Rolf Haaker

Senior Physician at the Orthopedics, Rheumatology, and Traumatology Clinic at St. Vincenz Hospital in Brakel

Frame orthosis

- + Effective guidance and stabilization according to the four-point principle, as well as flexibility during movement



Med. compression stocking

- + Effective post-operative and post-traumatic edema control using the VenoTrain soft S

BENEFITS OF TREATMENT

✓ **Improving rehabilitation**

✓ **Reducing existing edema**

✓ **Promoting healing**

THE PERFECT COMBINATION

SecuTec® Genu Flex

ENSURING INCREASED DYNAMICS DURING TREATMENT



- + Flexible during movement, thanks to unique Flex-Modules
- + Individual fit, adjustable over the course of treatment
- + Reliable protection during everyday activities

VenoTrain® soft S

THE ALL-ROUNDER FOR EASY DONNING OF THE PRODUCT



- + Anatomically contoured Easy Heel makes donning of the product easier
- + High level of wearing comfort around the foot, thanks to reduced knitted fabric in the instep area
- + Breathable, robust knitted fabric
- + Available as a knee-high or thigh-high in two compression classes as well as with an open or closed toe

For additional details about indications, contraindications, and handling of the required product, please consult the relevant instructions for use!



A complete overview of the orthoses portfolio can be found in the Orthopedic Product Catalog (item no. 190 0106 0052 102).





V6+ STANDARD SIZE CHART

		CIRCUM-FERENCE	I		II		III		IV		V		VI	
			NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
Length ℓ_T	85–108	c_T	67–89	71–91	73–111	75–115	76–118	79–119	82–127	84–129	86–136	91–138	96–140	102–145
Length ℓ_H	76–98	c_H	85–102		87–116		92–121		98–132		104–140		110–146	
Type	AT	c_G	45–53	51–59	48–56	54–62	51–59	57–66	54–63	60–69	57–66	63–73	60–70	66–76
Length ℓ_G		c_F	39–45	43–50	42–49	46–53	45–52	49–57	48–56	52–60	51–59	55–64	54–63	58–67
short	62–70	c_E	32–37	35–40	34–39	37–42	36–41	39–44	38–44	41–47	40–46	43–49	42–48	45–51
long	70–80	c_D	28–32	30–34	30–34	32–37	32–37	34–39	34–39	36–41	36–41	38–43	38–43	40–46
Type	AG AT	c_C	29–33	31–35	31–35	33–38	33–38	35–40	35–40	37–42	37–42	39–44	39–44	41–47
Length ℓ_D		c_{B1}	23–27	24–28	25–29	26–30	27–31	28–32	29–34	30–35	31–36	32–37	33–38	34–39
short	35–40	c_B	19–20,5		20,5–22		22–24		24–26		26–28,5		28,5–31	
long	40–44	c_Y	25–29		26–31		28–33		30–35		32–38		35–41	
Type	AD	c_A	19–22		20–23		22–25		24–27		26–30		28–33	

Sizing system information given in cm

VenoTrain® impuls

THE POWERFUL STOCKING

The VenoTrain impuls is a great accessory for patients suffering from pronounced venous insufficiency or soft connective tissue, or for those who prefer robust support for their active lifestyles. The stockings produce powerful deep action by exerting high working pressure with even the slightest muscle movements. Their microfiber material ensures strong support while maintaining wearing comfort.

Ideal for advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers and for early-stage lipedema and lymphedema.*

STANDARD COLORS



FASHION COLORS



FEATURES



- Particularly deep action and relief through high working pressure with even the slightest muscle movements
- High wearing comfort thanks to breathable microfiber material that is gentle on the skin
- Reduces constriction and creasing, especially in cases of soft connective tissue, thanks to short-stretch elastic knitted fabric
- Spacious Easy heel for easy donning and doffing
- Exceptional fit, thanks to highly differentiated V6+ sizing system (I–VI)

DETAILS

SIZING SYSTEM V6+ (I–VI)

COMPRESSION CLASSES

- Ccl 2 (23–32 mmHg)
- Ccl 3 (34–46 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 2: 80 % polyamide, 20 % elastane
- Ccl 3: 75 % polyamide, 25 % elastane

CARE INSTRUCTIONS



● STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG	AGHB	AT	ATH	ATX
Top band / Waist band**	● ▲			● ▲	▲	▲
Silicone dots	● ▲	● ▲				
Garter belt			● ▲			

* For a complete list of all indications, please see the VenoTrain instructions for use.

** For AT as standard

SELECTION


VenoTrain® SELECTION PORTFOLIO



Product guide

The VenoTrain Selection portfolio is an **exclusive collection of knee-high stocking solutions** that, in addition to classic compression stockings, also meet more individual demands. With their **exceptional materials** and an attractive design, the Selection products not only provide inconspicuous compression therapy without looking like compression garments, but also combine style with medical effectiveness.

This portfolio is versatile and suitable for different medical needs. It can be used for a wide range of venous disorders as well as other indications where compression therapy helps.



Choosing the right Selection product

1. indication and diagnosis provided by a physician
(product specified on prescription)
2. Compression class
(product specified on prescription)
3. Patient's individual requirements concerning material, functionality, and appearance

Color overview VenoTrain Selection portfolio

SIMPLE. CLEAR. CONSISTENT.

STANDARD COLORS		business	hero	cocoon
Black		●	●	●
Navy		●	●	●
Denim		●	●	●
Grey		●	●	●

SIZING SYSTEM

STANDARD SIZE CHART
PERFECT-FIT

DID YOU
KNOW?

This icon shows which products
are available with the
Perfect-Fit sizing system.



	CIRCUMFERENCE	S		M		L		XL	
		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
Length /D	cD	29–35	32–38	32–38	35–42	34–42	37–45	37–45	41–49
short 35–40	cC	31–37	34–40	34–40	37–44	37–45	39–47	40–48	43–51
long 40–46	cB1	25–29	27–31	27–32	29–34	29–35	31–37	32–38	34–40
xlong* 46–52	cB	20–22		22–25		25–28		28–31	
Shoe size	cY	29–33		31–35		34–40		37–44	
short 36–41	cA	19–24		22–27		25–31		28–35	
long 41–46									

* only available for VenoTrain hero
Sizing system information given in cm



VenoTrain® business

THE ELEGANT

While traveling or at work: the VenoTrain business combines style and elegance and fulfills all the requirements of a premium knee-high stocking.

With cotton used on the inside, our VenoTrain business is exceptionally breathable and skin-friendly, making it a perfect choice for business settings. The wide comfort band and reinforced sole provide additional wearing comfort.

FEATURES

- Gentle on the skin and comfortable thanks to the cotton material mix
- Reinforced sole for reliable foot support
- VenoTrain business retains its shape and appearance even after frequent wearing and washing thanks to a special fiber design with an anti-pilling effect

STANDARD COLORS

BLACK



NAVY



DENIM



GREY



DETAILS

SIZING SYSTEM PERFECT-FIT (S–XL)

COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed

MATERIAL CONTENT

- Ccl 1: 80 % polyamide, 14 % elastane, 6 % cotton
- Ccl 2: 71 % polyamide, 24 % elastane, 5 % cotton

CARE INSTRUCTIONS



● STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING



AD

Top band



Silicone dots



Available in shoe sizes 36–41 (foot short) and 41–46 (foot long)



MERINO
WOOL MIX



VenoTrain® hero

THE CASUAL

The VenoTrain hero is your perfect companion for every occasion, especially in your leisure time. While it shows its strengths during demanding professional activities, it also guarantees maximum comfort for restful moments.

Thanks to its seamless toe, the padded toe and heel sections with merino wool and a crease-free instep area, it offers a pleasant and comfortable feel throughout the day.

FEATURES

- **Three different leg lengths** to ensure the perfect fit
- Durable and hard-wearing thanks to robust materials
- Merino wool and fluffy padding in the foot area for maximum wearing comfort and a pleasant foot climate
- Also suitable for sensitive feet thanks to seamless toe
- Crease-free fit in the instep area thanks to special bending zone
- Cut-out at the Achilles tendon for targeted pressure relief



Seamless toe



Three different leg lengths

STANDARD COLORS

BLACK



NAVY



DENIM



GREY



DETAILS

SIZING SYSTEM PERFECT-FIT (S – XL)

COMPRESSION CLASSES

- Ccl 2 (23 – 32 mmHg)

TOE closed

MATERIAL CONTENT

- 84 % polyamide, 14 % elastane, 2 % cotton (Merino wool)

CARE INSTRUCTIONS



• STANDARD SIZES

TOP BAND EDGING



AD

Top band



Available in shoe sizes 36 – 41 (foot short) and 41 – 46 (foot long)





CLINICALLY
TESTED SKIN-
FRIENDLINESS



VenoTrain® cocoon

THE CARING

This medical compression stocking provides both compression and skin care, offering special protection and care for dry, sensitive skin.¹

High-quality care substances are integrated into the stocking to support the skin's natural barrier and make dry skin more supple.

STANDARD COLORS

BLACK



NAVY



DENIM



GREY



FEATURES

- Pleasant wearing comfort and optimum skin care thanks to innovative blend of natural fibers with integrated skin care substances
- Clinically proven to maintain skin hydration¹
- Maintains the skin condition during compression therapy by effectively combating itching and scaling

DETAILS



SIZING SYSTEM PERFECT-FIT (S – XL)

COMPRESSION CLASS

- Ccl 2 (23 – 32 mmHg)

TOE closed

MATERIAL CONTENT

- 52 % polyamide, 26 % elastane, 11 % cotton, 11 % cellulose

CARE INSTRUCTIONS



STANDARD SIZES

**TOP BAND
EDGING**




AD

Top band



Available in shoe sizes 36 – 41 (foot short) and 41 – 46 (foot long)



¹Westphal, T., et al. (2019) Medical compression stockings on the skin moisture in patients with chronic venous disease. Vasa doi: 10.1024/0301-1526/a000812

SPECIAL

VenoTrain® SPECIAL PORTFOLIO



Product guide

The VenoTrain Special portfolio meets specific requirements of patients with advanced vein problems, including venous leg ulcers (CEAP stages C3 – C6)*. With its compression garments, the portfolio also helps treat concomitant diseases, such as mild to moderate peripheral arterial disease, or PAD for short.



VenoTrain® angioflow

For compression therapy despite mild to moderate PAD¹.

VenoTrain® ulcertec

For the differentiated treatment of venous leg ulcers.

* For a complete list of all indications, please see the VenoTrain instructions for use.

¹ Stücker, M., et. al (2020). Sicherheit eines Kompressionsstrumpfes für Patienten mit chronischer venöser Insuffizienz (CVI) und peripherer arterieller Verschlusskrankheit (pAVK). JDDG, 18 doi: 10.1111/ddg.14042_g



SIZING SYSTEM

The VenoTrain ulcer tec is measured according to the **X5+ sizing system**.

NOTE:

Based on its special product features, the VenoTrain angioflow is **only** available as a custom-made product.

	CIRCUM-FERENCE	XS		S		M		L		XL	
		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
Length cG short 61–70 long 70–80 Type AG	cG	45–57	52–66	49–62	55–68	53–67	58–74	58–73	65–82	63–80	69–90
Length cD short 35–40 long 40–45 Type AD	cC	24–29	28–34	27–32	31–38	29–36	35–42	33–40	38–47	36–44	43–52
	cB	20–21,5		21,5–23		23–25		25–28		28–32	

All sizes in centimeters (cm)



VenoTrain® angioflow

THE SPECIALIST COMPRESSION PRODUCT FOR USE DESPITE PAD¹

The VenoTrain angioflow is a medical compression stocking designed with product features that allow it to be safely¹ used for compression therapy despite mild to moderate peripheral arterial disease (PAD) (ABI < 0.9 and > 0.5).

The specifically developed short-stretch elastic contoured knit ensures a high working pressure despite a low resting pressure. The deep action of the VenoTrain angioflow supports the venous pumping function and alleviates the symptoms of a venous disorder, thus having a positive effect on quality of life.

FEATURES

- Specifically developed contour knitting technique for individualized care and the perfect fit
- High working pressure during movement and low resting pressure in Ccl 1, thanks to short-stretch knitted fabric
- Low-compression comfort foot prevents pressure points based on the reinforced instep area and padded sole as well as the seam at the top above the toes
- Spacious Easy heel for easier donning and doffing
- Prevents constriction and offers a secure hold thanks to extra-wide silicone top band
- Skin-friendly microfiber material

▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING

Silicone dots



AD



COLORS

CARAMEL



DETAILS

SIZING SYSTEM Made-to-measure

COMPRESSION CLASS

- Ccl 1 (18–21 mmHg)

TOE closed

MATERIAL CONTENT

- 79 % polyamide, 21 % elastane

CARE INSTRUCTIONS



SUCCESSFUL TREATMENT BASED ON EXACT MEASUREMENTS

BODYTRONIC® 610



For successful treatment, the VenoTrain angioflow must fit like a second skin. It is therefore only available as a custom-made product. We recommend exact measurement using Bauerfeind's digital BODYTRONIC 610 measurement system. Only this measurement method allows precise custom-made "contour knitting". This is necessary for the manufacture of the special short-stretch knitted fabric.

¹ Stücker, M., et. al (2020). Sicherheit eines Kompressionsstrumpfes für Patienten mit chronischer venöser Insuffizienz (CVI) und peripherer arterieller Verschlusskrankheit (pAVK).

JDDG, 18 doi: 10.1111/ddg.14042_g

² Ankle-brachial index



C0

C1

C2

C3

C4

C5

C6

VenoTrain® angioflow

THE SPECIALIST COMPRESSION PRODUCT FOR USE DESPITE PAD¹

Medical compression stocking that can safely be used for compression therapy despite mild to moderate peripheral arterial disease (PAD)¹.

Based on their material properties, many compression stockings are unsuitable for compression therapy if patients also suffer from mild to moderate PAD. Low resting pressure tends to produce low working pressure, which is also called contact pressure, during movement because of the material's high elasticity. If a higher working pressure is generated, most circular knit compression stockings also have a higher resting pressure, meaning a higher compression class, based on their long-stretch characteristics.

Another characteristic of compression stockings is that they are able to withstand an increase in circumference, during movement for example. This is referred to as 'stiffness'. Underneath very stretchy materials, there is only a small increase in pressure (lower stiffness), while this pressure increase can be very high underneath less stretchy materials (higher stiffness).²

Compression therapy guidelines specify that a higher degree of stiffness is recommended for compression stockings used to treat conditions like severe chronic venous insufficiency (CVI) or peripheral artery disease.

STUDIES PROVE^{3,4}

With advancing age, the likelihood of suffering from chronic venous insufficiency (CVI) as well as peripheral arterial disease (PAD) increases.



¹ Stücker, M., et. al (2020). Sicherheit eines Kompressionsstrumpfes für Patienten mit chronischer venöser Insuffizienz (CVI) und peripherer arterieller Verschlusskrankheit (pAVK). JDDG, 18 doi: 10.1111/ddg.14042_g

² Medical Compression Therapy of the Extremities with Medical Compression Stockings (MCS), Phlebological Compression Bandages (PCB), and Medical Adaptive Compression Systems (MAC) – AWMF registration number: 037/005. Version 2023

³ Rabe, E. et al. (2003) Bonner Venenstudie der DGP, Phlebologie

⁴ Criqui, M. H., et al. (1985); The prevalence of peripheral arterial disease in a defined population. Circulation.

PROVEN BY RESEARCH¹

VenoTrain angioflow in patients with combined CVI / PAD after 14 days of wearing the product:

- + Constant arterial circulation
- + No negative skin changes (e.g. constriction, pressure points, skin damage)
- + Patients appreciate the product because of its easy handling and great fit
- + Patients rate the VenoTrain angioflow as very good to good
- + Significant alleviation of CVI problems (e.g. reduction in pain as well as in feelings of heaviness and swelling of the legs)

Product features ensuring the perfect care:

Low resting pressure Ccl 1 (18–21 mmHg)

- Suitable for the management of mild to moderate PAD (ABI < 0.9 and > 0.5), as the low resting pressure does not reduce arterial perfusion
- Noticeably reduced feeling of pressure during periods of inactivity, yet offering unmitigated medical effectiveness
- Easy donning and doffing

Short-stretch elastic knitted fabric: less elasticity means more resistance during movement

- Optimized massage effect and deep-acting support of the venous pumping function through high working pressure during movement
- Alleviation of the symptoms of venous disease, positively impacting quality of life

Available only as a custom-made solution:

- The perfect fit, thanks to the special contour knitting technique
- Stockings like a second skin – adapted to the patient's needs in cases of concomitant mild to moderate PAD



RHOMBOID
KNIT



VenoTrain® ulcertec

THE SPECIALIST COMPRESSION PRODUCT FOR VENOUS LEG ULCERS



For videos and more details about the VenoTrain ulcertec, please visit our online shop
experts.bauerfeind.com



The combination of an overstocking with a patented rhomboid knitted fabric specially developed by Bauerfeind and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.


The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

COLORS

BEIGE
OVERSTOCKING



WHITE
LINER



FEATURES

- Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort¹
- Differentiated treatment based on two designs within one compression class

39 – moderate

- Compression class 3 (34–46 mmHg)
- Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking

46 – strong

- Compression class 3 (34–46 mmHg)
- Liner with moderate compression providing a compression pressure in the upper range of compression class 3 in combination with the overstocking

DETAILS

SIZING SYSTEM X5+ (XS–XL)

COMPRESSION CLASSES

- 39–moderate: 34–46 mmHg
Liner with gentle compression
- 46–strong: 34–46 mmHg
Liner with moderate compression

TOE

- open (overstocking)
- closed (liner)

MATERIAL CONTENT

LINER AG²

- 39–moderate: 68 % polyamide, 32 % elastane

LINER AD

- 39–moderate: 87 % polyamide, 13 % elastane
- 46–strong: 69 % polyamide, 31 % elastane

OVERSTOCKING (39 – moderate / 46 – strong)

- Standard: 79 % polyamide, 21 % elastane
- Made-to-measure: 68 % polyamide, 32 % elastane

CARE INSTRUCTIONS

- Overstocking:
- Liner:

● STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING



AD



AG

Top band / Comfort band



Silicone sensitive



¹ Jünger, M., et al. (2004) Efficacy and tolerability of an ulcer compression stocking for therapy of chronic venous ulcer compared with a below-knee compression bandage: results from a prospective, randomized, multicentre trial.
Current medical research and opinion,
doi: 10.1185/030079904X4086

² The AG liner additionally has a silicone top band on the upper edge of the stocking. It consists of polyamide, elastane, and silicone.

BASIC

VenoTrain®

BASIC

PORTFOLIO

Product guide

The VenoTrain Basic portfolio provides cost-effective basic treatment of patients with chronic venous disorders.



SIZING SYSTEM

STANDARD SIZE CHART
PERFECT-FIT*

		CIRCUMFERENCE		S		M		L		XL	
				NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
		cT		up to 85		up to 90		up to 100		up to 120	
		cH		up to 95		up to 105		up to 120		up to 140	
Length IG		Silicone top band, wide		51–57	59–66	56–63	62–69	60–67	66–74	63–71	71–80
		Silicone top band, regular		45–51	52–59	49–56	55–62	53–60	58–66	55–63	62–71
short	61–70	Top band, AT, silicone top band		45–57	52–66	49–63	55–69	53–67	58–74	55–71	62–80
		cF		40–52	43–55	43–55	46–60	45–59	50–64	49–63	52–68
long	70–80	cE		33–40	36–44	35–43	39–49	38–48	42–53	42–53	47–58
Type	AG, AT	cD		29–35	32–38	32–38	35–42	34–42	37–45	37–45	41–49
		cC		31–37	34–40	34–40	37–44	37–45	39–47	40–48	43–51
Length ID		cB1		25–29	27–31	27–32	29–34	29–35	31–37	32–38	34–40
		cB		20–22		22–25		25–28		28–31	
short	35–40	cY		29–33		31–35		34–40		37–44	
		cA		19–24		22–27		25–31		28–35	
long	40–46										
Type	AD										

Sizing system information given in cm

For shoe sizes smaller than 36 or larger than 46 as well as in cases of particular anatomical features of foot, the use of custom-made products is recommended.

VenoTrain® pure

BASIC RANGE

The VenoTrain pure compression stocking supports venous activity and alleviates tired, heavy legs. Thanks to its skin-friendly and breathable knitted fabric, it provides a high level of wearing comfort.

FEATURES

- For men and women
- Available in the tried-and-tested Perfect Fit sizing system
- Hard-wearing material

COLORS

BEIGE



BLACK



DETAILS



SIZING SYSTEM PERFECT-FIT (S – XL)

COMPRESSION CLASSES

- Ccl 1 (18 – 21 mmHg)
- Ccl 2 (23 – 32 mmHg)

TOE closed, open

MATERIAL CONTENT

- Ccl 1: 77 % polyamide, 23 % elastane
- Ccl 2: 74 % polyamide, 26 % elastane

CARE INSTRUCTIONS



● STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

TOP BAND EDGING	AD	AG	AT
Top band / Comfort band	● ▲		● ▲
Silicone dots	● ▲	● ▲	

ACTIVE

VenoTrain[®] ACTIVE PORTFOLIO

The VenoTrain Active portfolio provides gentle support at work, during leisure time, and exercise. The VenoTrain act support stockings stand out with their revitalizing wearing sensation.


The Compression Socks are the perfect accessory for all sorts of outdoor and leisure activities.




SIZING SYSTEM

VenoTrain® act – SHEER ELEGANCE

AG, AT, ATU



AT CAL



AD

SIZE	0	1	2	3
SHOE SIZE	37–39	39–41	41–43	43–45

VenoTrain® act – COTTON AND MICROFIBER

AD

COTTON (for him)


SIZE	1	2	3	4
SHOE SIZE	39–41	41–43	43–45	45–47

AD

MICROFIBER (for men and women)

SIZE	0	1	2	3	4
SHOE SIZE	37–39	39–41	41–43	43–45	45–47

COMPRESSION SOCKS



VenoTrain[®] act SHEER ELEGANCE













PROVIDING ELEGANT SUPPORT

In 70 or 140 denier, the VenoTrain act sheer elegance provides light to medium support in a modern translucent design that is available in a wide range of colors and styles.

FEATURES

- For her
- Elegant and transparent
- Available in seven shades
- Available as a knee-high or thigh-high stocking with a silicone lace band and as pantyhose
- Also available as pantyhose for pregnant women with a special pantyhose top that expands as the mother's bump grows
- Calibrato pantyhose (AT CAL) for particularly feminine curves

STANDARD SIZES

							
			AD	AG	AT	ATU	AT CAL
DEN			70 / 140	70 / 140	70 / 140	140	70 / 140
COLORS		SAND	●	●	●		
		CARAMEL	●	●	●	●	●
		AMARETTO	●	●	●		
		ESPRESSO	●	●	●		
		BLACK	●	●	●	●	●
		GRAPHITE	●	●	●		
		NAVY	●	●	●		

VenoTrain act is not a medical device.

* We are happy to share example specifications, additional material contents with you on request.

No alterations, repairs, or additional services are offered for the VenoTrain act.
It cannot be returned in exchange for a credit note or another item.

experts.bauerfeind.com/venotrain-act-elegance

DETAILS

TOE closed

MATERIAL CONTENT*

AG (70 DEN) 81 % polyamide, 19 % elastane

AG (140 DEN) 73 % polyamide, 27 % elastane

CARE INSTRUCTIONS



SHEER ELEGANCE

VenoTrain[®] act COTTON

VenoTrain[®] act MICROFIBER

PROVIDING GENTLE SUPPORT

VenoTrain act cotton, available as a men's stocking, and VenoTrain act microfiber, available for men and women, offer a high degree of support in 280 denier and are reliable companions when traveling, at work or during leisure activities.

FEATURES






COTTON

- For him
- Particularly comfortable to wear, thanks to cotton
- With distinctive ribbed look
- Opaque
- Knee-high stocking with a high level of support (280 DEN)
- Available in four colors

MICROFIBER

- For men and women
- Particularly skin-friendly and soft, thanks to high microfiber content
- Fits snugly around the leg, while allowing the skin to breathe
- Knee-high stocking with a high level of support (280 DEN)
- Available in four colors

STANDARD SIZES

			 AD
DEN			280
COLORS		BEIGE	●
		BLACK	●
		GRAPHITE	●
		NAVY	●

DETAILS

TOE closed

MATERIAL CONTENT

- COTTON: 27 % polyamide, 10 % elastane, 63 % cotton
- MICROFIBER: 94 % polyamide, 6 % elastane

CARE INSTRUCTIONS



COTTON



MICROFIBER



COMPRESSION SOCK TRAINING

COMPRESSION SOCK PERFORMANCE

PROVIDING SUPPORT DURING EXERCISE

The **Compression Socks Training** and **Performance** were developed on the basis of extensive expertise in biomechanics and compression therapy.

They offer noticeable support during continuous exertion, prevent tiring muscle vibrations and help the muscles to regenerate more quickly.

FEATURES

- Muscle Toning Zone supports the calf muscles
- Shock absorption zones protect the forefoot and heel areas
- Mild compression to stimulate the foot muscles
- The taping zone helps prevent injuries and is designed to improve proprioception
- Especially suited for sports with stopping and jumping movements



TAPING ZONE

PERFORMANCE

- Designed to improve muscle efficiency
- Protection zone on the Achilles tendon
- Breathable microfiber
- Especially suitable for endurance sports not requiring many changes of direction

DETAILS

TOE closed

MATERIAL CONTENT

TRAINING 87 % polyamide, 13 % elastane

PERFORMANCE 75 % polyamide, 25 % elastane

CARE INSTRUCTIONS

TRAINING

PERFORMANCE

TRAINING AVAILABLE IN



Charcoal / Rivera

Charcoal / Polar

Silver / Rivera

Silver / Polar

PERFORMANCE AVAILABLE IN



Rivera

Black

White

Wild
Berry

French
Lavender

NEW
SPECIAL
EDITION

STANDARD SIZES

TOP BAND EDGING

Top band /
Comfort band

TRAINING

PERFORMANCE



AD



Compression Sock Training and Performance are no medical devices. No alterations, repairs, or additional services are offered for the Compression Socks. It cannot be returned in exchange for a credit note or another item.

experts.bauerfeind.com/compression-sock-training
experts.bauerfeind.com/compression-sock-performance

OEKO-TEX®
CONFIDENCE IN TEXTILES
STANDARD 100
598-0153 HOFENSTEIN HTTL
Tested for harmful substances.
www.oeko-tex.com/standard100

ADDITIONAL PRODUCTS THE PERFECT ADDITION FOR EVERYDAY USE

The functional additions for medical compression therapy with VenoTrain by Bauerfeind enhance the ease of use of compression stockings. They provide care for both the skin and the compression stockings, thus supporting and facilitating treatment.





**TO BE APPLIED IN THE
MORNING AND EVENING
WHEN USING MEDICAL
COMPRESSION STOCKINGS**

BAUERFEIND CARE FOAM*

- Contains panthenol to protect and soothe particularly stressed and sensitive skin
- Before putting on the stockings and in the evening after taking them off



LINGERIE BAG

To protect the compression stockings during every machine wash



DETERGENT*

Cleans gently, ensuring careful product treatment and providing lasting protection



**FOR USE IN THE
EVENING AFTER
DOFFING OF
YOUR MEDICAL
COMPRESSION STOCKINGS**

VENOTRAIN® LOTION*

- Developed in collaboration with Sebapharma
- Contains chamomile extract and is specially designed for people with sensitive skin
- For use in the evening after taking off their stockings



**BAUERFEIND
SPECIAL GLOVES**

Help with putting the compression stockings on correctly

VenoTrain® glider

COMPACT DONNING AID

A video with donning instructions can be found at:
bfzr.de/glider-plus



The VenoTrain glider is a practical aid for putting on medical compression stockings.

The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe.

Saving space and time, it is the ideal companion when you're out and about.

FEATURES

- Simple handling of the compression stockings for easier donning
- Also ideal for traveling thanks to its extremely lightweight design
- Safe application thanks to sturdy and tear-resistant material
- Product is available on prescription



DONNING OF THE STOCKINGS: STEP 1



DONNING OF THE STOCKINGS: STEP 2



DONNING OF THE STOCKINGS: STEP 3



DONNING OF THE STOCKINGS: STEP 4



DONNING OF THE STOCKINGS: STEP 5

VenoTrain® glider plus

ROBUST DONNING AID

The VenoTrain glider plus is a robust donning aid with an ergonomic design for all types of medical compression stockings.

Folding handles make it easy to put on the stockings, avoiding the need to bend down.

FEATURES

- Simple handling and effort-saving, ergonomic design of the compression stockings for easier donning
- Safe application using the lightweight metal frame and the two folding handles on the sides
- Compact design for space-saving storage
- Product is available on prescription



DONNING OF THE STOCKINGS: STEP 1



DONNING OF THE STOCKINGS: STEP 2



DONNING OF THE STOCKINGS: STEP 3



DONNING OF THE STOCKINGS: STEP 4



DONNING OF THE STOCKINGS: STEP 5



DONNING OF THE STOCKINGS: STEP 6

BAUERFEIND MEASUREMENT TECHNOLOGY

QUALITY THAT CAN BE MEASURED –
AND MADE-TO-MEASURE.



BODYTRONIC® 410 FLEXIBLE 3D SCANNER FOR THE LOWER EXTREMITIES

- Contact-free measurement for compression stockings, with no need for special measurement clothing
- 3D image of the measurement area up to the waist
- Quick and easy to set up and take down



BODYTRONIC® 610 FULL-BODY SCANNER INCLUDING PRODUCT SUPPLY FOR ALL EXTREMITIES

- Contact-free measurements for medical compression stockings, knee and foot supports, knee orthoses as well as for back, armband and shoulder products
- Optional full-body scan, export of 3D data
- Standing space of only 1.60 m²

BAUERFEIND ACADEMY

COMPETENT. GOAL-ORIENTED. SUCCESSFUL.

Would you like customers to clearly favor your medical retail store and reward your team's expertise?

We will be happy to support you!

We offer you customized training with strong practical relevance. These training opportunities can also be conducted online.

Discuss the options with your sales manager.

BENEFITS AT A GLANCE

+ TRAINING WITH STRONG PRACTICAL RELEVANCE

The training sessions offered by our Academy will always include medical information to improve understanding and ensure increased technical expertise.

Specifically trained coaches, such as orthotists, physical therapists, and lymph therapists, will provide professional guidance through the different subjects, while always considering the relevant target group.

+ IN-HOUSE TRAINING SESSIONS ACCORDING TO YOUR SPECIFICATIONS

Are you interested in personalized training that is tailor-made to your specific requirements?

You can book your training session according to your specifications. You choose the subject, the location and the time. Just get in touch!

+ ONLINE TRAINING

Would you like to learn in a time-saving and cost-efficient way? You can book tailored online training and specify the contents and the subject according to your needs.

+ E-LEARNINGS

Explore our e-learning modules from the Bauerfeind Academy – learn at your own speed, wherever you are.

bfzr.de/ph-elearning



For further information and registration forms, please contact your local Sales Representative.

SERVICE GUARANTEES

INTERNATIONAL CUSTOMER SERVICE



100 % AVAILABILITY

By telephone within our service times
(Mo. – Th. 7.00 am – 05.00 pm,
Fr. 7.00 am – 04.00 pm CET).

By fax, e-mail and electronic ordering
systems 24 hours per day.



ELECTRONIC ORDERING SYSTEM

Make use of our offers for a fast,
faultless and effort reduced data
transfer: online shop, EDI service partner.



DISPATCH HANDLING

At your request we take care of the
full processing of your shipment
(incl. shipment quotations, additional
documents) and we can bundle your
order flexibly.



PROCESSING OF RETURNS

Within 5 days after receipt of goods.



ONE FACE TO THE CUSTOMER

You have one dedicated contact person
for all of your desires and concerns.



PROACTIVE INFORMATION

We keep you posted – you will receive
every important information about your
order and delivery status.



RETURN OF GOODS

Redemption of faultless goods within
30 days after delivery with 100 % refund
or replacement (T&C's apply).



BAUERFEIND- ONLINE SHOP




ONLINE SHOP SERVICE

When ordering via the Bauerfeind online shop, you benefit from reduced shipping costs and can take advantage of versatile functions such as the order status, availability display or patient files. To create a customer account, register at:

experts.bauerfeind.com/register

Order online now:
experts.bauerfeind.com



YOUR BENEFITS AT A GLANCE

- + 24 h online availability
- + Accessible anywhere using a web browser
- + Intuitive order process
- + Up-to-date order status
- + Free shipping for custom-made solutions
- + Smart search function
- + Digital order process
- + Direct connection to Bauerfeind BODYTRONIC measurement technology
- + Patient file to store measurement and treatment data
- + Support via Customer Service, chat, e-mail, contact form:
experts.bauerfeind.com/contact

EDGING AND TOP BANDS

PERFECT HOLD WITH A HIGH LEVEL OF COMFORT

EDGING

SECURE AND COMFORTABLE




Similar image

TOP BAND

For all knee-high compression stockings

The soft edging ensures wearing comfort and doesn't dig into the skin. Its neutral look works for both women and men.

Our edging and top bands make sure that compression stockings fit securely and don't slip. Their soft materials are breathable, gentle on the skin and very comfortable to wear. Learn more about the individual styles and different designs.



Similar image

COMFORT BANDS

For all pantyhose

This edging provides an excellent fit and a high level of comfort thanks to its higher elasticity and wide band. It sits comfortably without constricting.



WAIST BAND

For VenoTrain impuls pantyhose (standard)

Flat edging with an elastic band without silicone coating that offers exceptional wearing comfort.



GARTER BELT

The garter belt without silicone coating improves the hold of the stockings and offers plenty of freedom of movement. It is suitable for sensitive skin and can be worn on one or both sides.

SILICONE TOP BANDS

SECURE AND INDIVIDUALIZED



SILICONE MICRO DOTS

The highly elastic silicone top band with its small, flat silicone dots is very inconspicuous, even when worn under tight clothing, and leaves no marks on the skin. The skin-friendly material with its soft surface perfectly adapts to the contours of the leg.



LACE SENSITIVE

The highly elastic silicone top band has an ultra-thin, continuous coating that provides a perfect hold. It adapts perfectly to the leg contour, even with soft connective tissue.



SILICONE DOTS

Our classic silicone top band is mainly used for thigh-high stockings. Skin-friendly material, neutral design. For all those who like it simple.



SILICONE DOT PATTERN

The skin-friendly silicone top band with an elegant pattern.

STANDARD PRODUCT RANGE

CIRCULAR-KNIT PRODUCTS


DESIGN		SIZE	CCL			EDGING / TOP BAND								TOE	
			1	2	3	Top band / Comfort band	Waist band	Lace sensitive	Silicone micro dots	Silicone dots	Silicone dot pattern	Silicone sensitive	Garter belt	closed	open
VenoTrain															
discrétion	AD	I – IV	●*	●**		●								●	●
	AG	I – IV	●*	●**						●	●			●	●
	AT	I – IV	●*	●**		●								●	●
micro	AD	S – XL	●	●		●				●				●	●
	AG	S – XL	●	●				●	●	●	●			●	●
	AT	S – XL	●	●		●								●	●
	ATU	S – XL	●	●		●								●	●
micro Design Edition (Jive, Tango)	AD	S – XL	●	●		●								●	
	AG	S – XL	●	●				●	●	●	●			●	
	AT	S – XL	●	●		●								●	
soft	AD	S – XL	●	●		●				●				●	●
	AG	S – XL	●	●				●	●	●	●			●	●
	AGHB	S – XL	●	●									●	●	●
	AT	S – XL	●	●		●								●	●
soft S	AD	XS – XL	●	●		●				●				●	●
	AG	XS – XL	●	●				●	●	●	●			●	●
impuls	AD	I – VI		●	●	●				●				●	●
	AG	I – VI		●	●					●				●	●
	AGHB	I – VI		●	●								●	●	●
	AT	I – VI		●	●		●							●	●
business	AD	S – XL	●	●		●								●	
hero	AD	S – XL		●		●								●	
cocoon	AD	S – XL		●		●								●	
ulcertec 39 moderate	AD		34 – 46 mmHg			●								Liner	Over- stocking
	AG	XS – XL										●			
	AD/AG					●									
ulcertec 46 strong	AD	XS – XL	34 – 46 mmHg			●								Liner	Over- stocking
pure	AD	S – XL	●	●		●				●				●	●
	AG	S – XL	●	●						●				●	●
	AT	S – XL	●	●		●								●	●

* 10 – 15 mmHg ** 15 – 20 mmHg

VERSIONS



AD
Knee-high
stocking



AF
Mid-high
stocking



AG
Thigh-high
stocking



AGHB
Thigh-high
compression
stocking with
garter belt



AT
Pantyhose



ATU
Maternity
pantyhose



ATH
Gent's
pantyhose
without fly



ATX
Pantyhose
with com-
pressive
pantyhose
top



ATE
Single-leg
pantyhose



Toe
closed



Toe
open

MADE-TO-MEASURE PRODUCT RANGE

CIRCULAR-KNIT PRODUCTS

DESIGN		CCL			EDGING / TOP BAND								TOE			
		1	2	3	Top band / Comfort band	Lace sensitive	Silicone micro dots	Silicone dots	Silicone dot pattern	Silicone sensitive	Garter belt		closed	open	No foot	Stump
VenoTrain																
	micro															
	AD	▲	▲		▲			▲					▲	▲		
	AF	▲	▲			▲	▲	▲	▲				▲	▲		
	AG	▲	▲			▲	▲	▲	▲				▲	▲		
micro	AT	▲	▲		▲								▲	▲		
	ATU	▲	▲		▲								▲	▲		
	AD	▲	▲		▲								▲			
	AG	▲	▲			▲	▲	▲	▲				▲			
	AT	▲	▲		▲								▲			
soft	AD	▲	▲		▲			▲					▲	▲	▲	▲
	AF	▲	▲			▲	▲	▲	▲				▲	▲	▲	▲
	AG	▲	▲			▲	▲	▲	▲				▲	▲	▲	▲
	AGHB	▲	▲								▲		▲	▲	▲	▲
	AT	▲	▲		▲								▲	▲	▲	▲
	ATU	▲	▲		▲								▲	▲	▲	
	ATH	▲	▲		▲								▲	▲	▲	
	ATX	▲	▲		▲								▲	▲	▲	
	ATE	▲	▲		▲								▲	▲	▲	
	AD	▲	▲		▲			▲					▲	▲		
	AG	▲	▲			▲	▲	▲	▲				▲	▲		
	AGHB		▲	▲				▲					▲	▲		
	AT		▲	▲	▲								▲	▲		
impuls	AD		▲	▲	▲			▲					▲	▲		
	AG		▲	▲				▲					▲	▲		
	AGHB		▲	▲							▲		▲	▲		
	AT		▲	▲	▲								▲	▲		
	ATH		▲	▲	▲								▲	▲		
	ATX		▲	▲	▲								▲	▲		
business	AD	▲	▲		▲			▲					▲			
angioflow	AD	▲						▲					▲			
ulcertec 39 moderate	AD	34–46 mmHg			▲								Liner	Over- stocking		
	AG									▲						
ulcertec 46 strong	AD	34–46 mmHg			▲								Liner	Over- stocking		
pure	AD	▲	▲		▲			▲					▲	▲		
	AG	▲	▲					▲					▲	▲		
	AT	▲	▲		▲								▲	▲		



ADDITIONAL SERVICES ALTERATIONS | REPAIRS

INDIVIDUAL CUSTOMER SERVICE

ADDITIONAL SERVICES FOR NEW ORDERS

▲ MADE-TO-MEASURE SOLUTIONS

	micro	micro DESIGN EDITION	soft	soft S	impuls	business
Knitted toe area To be lengthened or shortened (3 – 10 cm)	▲		▲	▲	▲	
Toe Seam on the outside	▲	▲	▲	▲	▲	▲
Toe Wider, narrower, anatomical	▲	▲	▲	▲	▲	▲
Panty area Front/back – two seams	▲		▲*		▲	
ATH with longitudinal slit			▲		▲	
ATH with transverse opening			▲		▲	
AT with extra-wide gusset for the fuller figure	▲		▲		▲	
Sew in a stoma opening (only for repeat order)			▲			

* Please get in touch with the Bauerfeind Customer Service Team.

ALTERATIONS AND REPAIRS

● STANDARD SIZES

▲ MADE-TO-MEASURE SOLUTIONS

		micro	micro DESIGN EDITION	soft	soft S	impuls	business	hero	cocoon	ulcertec	pure
Top of stocking	Replace silicone top band	● ▲	● ▲	● ▲	● ▲	● ▲	▲			● ▲	● ▲
	Change top band to silicone top band	● ▲		● ▲	● ▲	● ▲	▲				
	Silicone top band narrower or wider	● ▲	● ▲	● ▲	● ▲	● ▲	▲			● ▲	● ▲
	Remove hip section and sew on silicone top band			● ▲		● ▲					
AT	Shorten pantyhose top (except ATX)	● ▲	● ▲	● ▲		▲					● ▲
	Lengthen pantyhose top			▲ *		▲					
	Make gusset narrower (excluding ATX)	▲ *		▲ *		▲					▲
	Make gusset wider (excluding ATX)	▲ *		▲ *		▲					▲ *
	Extra-wide gusset for the fuller figure (excluding ATX)	▲		▲		▲					
	Alter AT to AG with silicone top band	▲	▲	▲		▲					▲
AGHB	Replace hip strap			● ▲		● ▲					
	Shorten hip section			● ▲		● ▲					
	Lengthen / shorten hip strap			▲		▲					
Add a stoma opening or an open crotch afterwards (except ATX)				▲							
Shorten stocking	AD by up to 3 cm	● ▲	● ▲	● ▲	● ▲	● ▲				● ▲	● ▲
	AF by up to 4 cm	▲		▲							
	AG by up to 5 cm	● ▲	● ▲	● ▲	● ▲	● ▲				● ▲	● ▲
Redye from light to dark**		▲	▲	▲	▲	▲					
Toe alteration	Alter closed toe to open toe	▲ *		▲ *	▲ *	▲ *					▲ *
	Shorten open toe (from 3 cm)	▲		▲	▲	▲				▲	▲
Repairs, if possible	Hole*	● ▲		● ▲	● ▲	● ▲	● ▲	●	●	● ▲	
	Run*	● ▲		● ▲	● ▲	● ▲	● ▲	●	●	● ▲	
	Panty area seam	● ▲		● ▲		▲					
	Loose thread*	● ▲		● ▲	● ▲	● ▲	● ▲	●	●	● ▲	

*Please get in touch with the Bauerfeind Customer Service Team.

**Please do not try on stockings – re dyeing is no longer possible after contact with the skin.

Please refer to the current price list for information about designs, item numbers and prices for standard and made-to-measure solutions. If you have any questions about our products, please contact the Bauerfeind Customer Service Team.

MEASURING CORRECTLY FOR OPTIMIZED TREATMENT SUCCESS

Medical compression stockings can only be fully effective when they fit properly. This requires, for example, the exact measuring of all relevant body dimensions. Our order pads will make your daily work easier.

We are happy to send you order pads free of charge.

The different order sheets can also be downloaded at any time:

experts.bauerfeind.com/downloads



VENOTRAIN ORDER PAD

Item number

Order sheet for custom-made products on the front, order sheet for standard products on the back	271 2020 0000 002
--	-------------------

ORDER PAD FOR CUSTOM-MADE CIRCULAR KNIT PRODUCTS

Item number

VenoTrain soft arm	271 1220 3020 001
--------------------	-------------------

ADDITIONAL BAUERFEIND PRODUCT CATEGORIES



A healthy body needs healthy feet. **ORTHOPEDIC FOOT ORTHOSES** simultaneously relieve the feet, knees, hips, and back. They help the foot to function naturally and can therefore positively influence the body's biomechanics.

As different feet need different types of support, Bauerfeind offers more than 400 different types of foot orthosis.

**MEDICAL AIDS
HAVE BEEN
PROVEN TO BE
EFFECTIVE***



SUPPORTS AND ORTHOSES

by Bauerfeind help after injuries as well as in cases of acute and chronic orthopedic conditions.

They improve the feeling of stability and provide confidence during everyday activities and exercise. The supports and orthoses are made from lightweight and breathable materials to perfectly underpin therapeutic exercise.

They are easy to handle during everyday use, thanks to innovative solutions such as quick-fit fasteners and integrated donning aids, and will win you over with premium quality.



COMPRESSION STOCKINGS FOR EDEMA TREATMENT

Every day is more pleasant when you're feeling good. That is why the flat knit, custom-made VenoTrain® curaflow compression stockings were developed.

This versatile and colorful product perfectly adapts to the individual shape of the lymphedema or lipedema patient's legs, arms, hands, or thorax.

VenoTrain® curaflow alleviates pain caused by tension and provides noticeable relief so patients can live their everyday lives feeling carefree, confident, and attractive.

*** Supports / orthoses, compression stockings and foot orthoses help patients get back on their feet quickly after an injury or in case of disorders.**

(Source: representative patient survey conducted by the Allensbach Institute for Public Opinion Research on behalf of eurocom e.V., 2023, n=1,300 users of supports/orthoses, compression stockings and foot orthoses)

BAUERFEIND AG TERMS AND CONDITIONS OF SALE AND DELIVERY (GTC)

Prices

The prices quoted are net prices in euros ex works and do not include statutory VAT.

Shipping

Shipping is at the buyer's cost and risk even if the buyer does not bear the transportation cost. Disruptions of operations affecting us or any of our suppliers, lack of suitable personnel, strike, lockout, and force majeure are causes for extending delivery and relieve us of any delivery obligation for the time of the obstruction.

Shipping costs

We calculate our shipping charges as a flat rate of EUR 7.00 for a net total order value of up to EUR 50.00, and EUR 4.00 for a net total order value of up to EUR 200.00. No shipping costs applied, when the order contains custom-made goods. Special shipping requirements will be charged at cost.

Returns

Returns outside the warranty terms are only accepted after prior advice and approval from us, provided the original delivery documents are submitted; it is supposed that the returned goods are in good order. We reserve the right of rejecting the acceptance of any returned consignment. Except where covered by warranty, we charge for processing returns. For returns within one month of the date of delivery, this charge is EUR 10.00 for each return consignment. For goods returned one month or more after the date of delivery, a handling fee of 20 percent of the net order value, or a maximum of EUR 25.00 per product and unit, is charged; for returns six to twelve months after delivery date, the fee is 40 percent of the net order value, or a maximum fee of EUR 50.00 for each product and unit.

Goods returned outside the warranty or 12 months or more after the date of delivery will not be accepted back. Custom made goods are non-returnable.

Payments

If not otherwise agreed individually in advance and in writing, Payments are due immediately following delivery. We offer a 4 percent discount for payments made by direct debit and a 3 percent discount for payments made within 10 days of the invoice date. Otherwise, the period allowed for payment is 30 days net. The buyer shall be considered in default once 30 days have passed since the invoice date. From this point on, any further deliveries will only be made against advance payment. From the time at which payment is delayed, we will charge default interest of eight percent above the applicable base interest rate. Payments by check, bank draft, etc. are deemed to have been made when the amount has finally been credited to our account.

Retention of title

We retain the title to all goods delivered until the buyer has paid all amounts outstanding under our business relations with him. The retention of title shall also apply proportionately to goods which have been processed or otherwise modified. If the goods are sold, all third-party claims are deemed to have been assigned to us without provision up to the total amount owing to us.

Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment we can demand that the buyer discloses the assigned claims and the names of the debtors, provides all information required for collection of the debt, delivers all required documents

and informs the debtors (third party) of the assignment. The buyer shall not be entitled to sell, pledge or assign the goods as security. Pledges, confiscation, or any other dispositions taken by a third party shall be promptly communicated to us by the buyer, and all information and documents necessary for us to defend our rights shall be provided. We undertake to release, on the buyer's demand, any security which exceeds a secured claim that has not been satisfied by more than 20 percent

Complaints

Warranty claims in respect of obvious defects shall be void if we do not receive immediate written notification of these claims in accordance with Section 377 of the German Commercial Code (HGB). Warranty claims shall also lapse one year after the delivery of the purchased goods. If the product is defective, we shall make a replacement delivery or repair it at our discretion. For reasons of hygiene, repairs will only be performed on products that have been washed or are completely hygienic. Upon receipt of unwashed goods, we reserve the right to claim a suitable sum for additional work. The customer may not demand a price reduction or withdraw from the contract provided that we fulfil our duties of subsequent performance, particularly with regard to the repair of defects or the delivery of goods without defect, unless our attempts at subsequent performance have failed. Defects which are caused by normal wear and tear, overuse, misuse, neglect of care, or noncompliance with the instructions for use are not covered by warranty.

Liability

Liability for damages is excluded unless there has been grossly negligent or willful conduct on our part or on the part of our representatives or agents. This limitation of liability shall not apply to damage due to injury to life, body or health, or the violation of a material contractual duty.

Place of jurisdiction

The place of performance for the delivery is our company headquarters. The place of jurisdiction for all disputes arising from the business relation is Gera, Germany.

Severability clause

The event that any individual provision of these terms and conditions of sale and delivery is or becomes invalid, will not affect the validity of the remaining provisions. The invalid provision shall be replaced with a valid provision which reflects the original economic intention of the parties as closely as possible.

Data protection clause / recording of conversations

For quality- and other purposes, we record calls, faxes, e-mails and letters between the buyer and our employees. All recordings serve legal provisions, which have to be met with medical products ordered. Recordings will be deleted automatically in line with legal provisions. Non of the recorded data will be disclosed to third parties, if not required by decision of court or on request of local tax or healthcare authorities. The buyer hereby consents to recordings such.

Please see www.bauerfeind.de/ds for further information on the latest version of Bauerfeind AG's General Data Protection Policy.

Last update: February 2022

For additional details about indications, contraindications, and handling of the required product, please consult the relevant instructions for use.

DISCLAIMER:

All information given in this catalog is up to date at the day of first publishing. The catalog allows – depending on the information given – an initial overview and suggests some, but not all options of medical diagnosis and treatments to experts only. This catalog does not address patients. Although scientifically based, the catalog by no means claims to be complete or show all existing doctrinal. The entire content of this catalog – also in extracts, no matter how published, and especially concerning diagnostics and forms of treatment (in text, pictures,

graphics, images or other material) has to be understood as non-binding support for the therapy to be selected finally by medical doctor, physician or other medical professionals. The catalog's intension is not to offer the full range of all medical aspects and latest scientific discussion that should be mandatorily considered before taking actions. We strongly recommend individual medical advice. If in doubt, please seek additional opinion of medical experts, who we would be delighted to provide.

MOTION IS LIFE

WITH COMPRESSION BY BAUERFEIND

Do you have any questions, suggestions or requests?
Please contact us:

BAUERFEIND AG

Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany
P +49 (0) 36628 66-10 00
F +49 (0) 36628 66-19 99
E info@bauerfeind.com

AUSTRIA

Bauerfeind Ges.m.b.H.
Hainburger Straße 33
1030 Wien
P +43 (0) 800 4430-130
F +43 (0) 800 4430-131
E info@bauerfeind.at

BENELUX

Bauerfeind Benelux B.V.
Waarderveldweg 1
2031 BK Haarlem

THE NETHERLANDS

P +31 (0) 23 531-9427
F +31 (0) 23 532-1970
E info@bauerfeind.nl

BELGIUM

P +32 (0) 2 527-4060
F +32 (0) 2 792-5345
E info@bauerfeind.be

BOSNIA AND HERZEGOVINA

Bauerfeind d.o.o.
Meše Selmovića 19
71000 Sarajevo
P +387 (0) 33 710-100
F +387 (0) 33 619-422
E info@bauerfeind.ba

CROATIA

Bauerfeind d.o.o.
Goleška 20
10020 Zagreb
P +385 (0) 1 6542-855
F +385 (0) 1 6542-860
E info@bauerfeind.hr

FRANCE

Bauerfeind France S.A.R.L.
13 rue de la Perdrix
93290 Tremblay-en-France
P +33 (0) 1 4863-2896
F +33 (0) 1 4863-2963
E info@bauerfeind.fr

ITALY

Bauerfeind AG
Filiale Italiana
Via Cornaggia, 58
20092 Cinisello Balsamo MI
P +39 02 8977 6310
F +39 02 8977 5900
E info@bauerfeind.it

REPUBLIC OF NORTH MACEDONIA

Bauerfeind Dooel Skopje
50 Divizija 24 a
1000 Skopje
P +389 (0) 2 3179-002
F +389 (0) 2 3179-004
E info@bauerfeind.mk

SERBIA

Bauerfeind d.o.o.
102 Omladinskih brigada
11070 Novi Beograd
P +381 (0) 11 2287-050
F +381 (0) 11 2287-052
E info@bauerfeind.rs

SINGAPORE

Bauerfeind Singapore Pte Ltd.
Blk 41 Cambridge Road #01-21
Singapore 210041
P +65 6396-3497
F +65 6295-5062
E info@bauerfeind.com.sg

SLOVENIA

Bauerfeind d.o.o.
Dolenjska cesta 242 b
1000 Ljubljana
P +386 (0) 1 4272-941
F +386 (0) 1 4272-951
E info@bauerfeind.si

SPAIN

Bauerfeind Ibérica, S.A.
Calle Albacete, 23,
Entresuelo C
46007 Valencia
P +34 96 385-6633
F +34 96 385-6699
E info@bauerfeind.es

SWITZERLAND

Bauerfeind AG
Vorderi Böde 5
5452 Oberrohrdorf
P +41 (0) 56 485-8242
F +41 (0) 56 485-8259
E info@bauerfeind.ch

UNITED ARAB EMIRATES

Bauerfeind Middle East FZ LLC
Dubai Healthcare City
Building 40, P.O. Box 505116
Dubai
P +971 4 4335-684
F +971 4 4370-344
E info@bauerfeind.ae

USA

Bauerfeind USA, Inc.
1230 Peachtree St NE
75 14th St NE
Suite 2350
Atlanta, GA 30309
P +1 800 423-3405
P +1 404 201-7800
F +1 404 201-7839
E info@bauerfeindusa.com